Italiano

Scopo

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English

Intended use

The Free Light Chains APS Kit, when used in conjunction with Beckman-Coulter "APS-Array®" Analyzer, is intended for the determination of *Free Light Chains (FLC) kappa and lambda (separately) in unconcentrated urine*.

Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with absorbed antibodies specific to "hidden" sites of Light Chains of Immunoglobulins.

Contents

Kit Type				K
<i>Antiserum Reagents</i>				
<input type="checkbox"/> Free Light Chains kappa	<i>R.APS.FRK</i>	green	1 vial – 4,5 ml	•
<input type="checkbox"/> Free Light Chains lambda	<i>R.APS.FRL</i>	yellow	1 vial – 4,5 ml	•
<i>Calibrators</i>				
<input type="checkbox"/> Free Light Chains kappa	<i>W.UPC.FRK</i>	<i>green</i>	1 vial – 1 ml	•
<input type="checkbox"/> Free Light Chains lambda	<i>W.UPC.FRL</i>	<i>yellow</i>	1 vial – 1 ml	•
<input type="checkbox"/> Prediluted Calibrator Set FLC Set of ready-to-use Calibrators	<i>W.SET.UPC.FRK.FRL</i>		5 vials – 1 ml 5 vials – 1 ml	-
<i>Controls</i>				
<input type="checkbox"/> Free Light Chains kappa Control	<i>Y.UPR.FRK</i>		1 vial – 1 ml	-
<input type="checkbox"/> Free Light Chains lambda Control	<i>Y.UPC.FRL</i>		1 vial – 1 ml	-
• Reagents included in the kit.				
▯ Correlated reagents not included in the kit. They are available separately.				

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

Information re Calibrators / Controls

Due to lack of provisions from the Scientific Committees our "Internal Standardisation" has been applied.

See the relevant leaflet enclosed.

Analysér Set up

Follow the manufacturer's instructions to program two "UDR" assays.

The assay parameters can be found on the back.

Processing Procedure

Calibration Curve

Use the Prediluted Calibrator Set or prepare various calibrator levels as shown in the table on the back.

Calibration Procedure

The calibration procedure must be run for each chemistry.

The calibration curve must be rerun at every cartridge replacement and whenever a new lot number is used

Put the concentrations of the calibration points on axis-x and the corresponding "Results" (Reagent Blank net) on axis-y. Results and curves are reported on the back.

Sample

Unconcentrated urine, centrifuged (3,000 rpm x 10'). Use 0,1% Sodium Azide as preservative.

Assay on samples

Samples must be run with FLC kappa and FLC lambda reagents.

To calculate sample concentration interpolate sample Rate result (Reagent Blank net) on the Calibration curve.

Quality Control

For each set of determinations it is recommended to include:

<input type="checkbox"/> Negative Control	PBS (or Saline) = Reagent Blank.
<input type="checkbox"/> Positive Control	<input type="checkbox"/> Calibrator dilution 1:16 - Conc. ≈ 1.25 mg/dl. The calibrator dilution is stable for a few hours. <input type="checkbox"/> As a Positive Control alternative use: Y.UPR.FRK and Y.UPR.FRL

Results Interpretation

Reference Range

Free Light Chains are virtually absent in the urine of a healthy subject.

Significant Values

Results higher than PBS results (reagent blank, negative control) are to be considered potentially positive. Double PBS amount can be considered as cut-off. 1 mg/dl is generally assumed as cut-off for screening purposes.

⚠️ Limitations – lack of sample/calibrator parallelism

The Free Light Chains in the sample could be "different" from those in calibrators and therefore they could react unproportionately; *lack of sample/calibrator parallelism*. The same occurs with the quantification of serum monoclonal immunoglobulins (serum Monoclonal Component – MC) by liquid phase ImmunoPrecipitation (turbidimetry, nephelometry).

⚠️ Antigen Excess

The application allows a precise assay of a calibrator concentration ten times higher than the highest set point of the calibration curve. However, Antigen Excess cannot be ruled out.

⚠️ Monoclonality

This test does not allow highlighting if the Free Light Chains in the sample are monoclonal or polyclonal.

⚠️ Advice on Reagent usage

100 tests can be run on one reagent vial. Reagent vial dead volume is 400 µl.

⚠️ Precautions and Hazards

Preservative: *Sodium Azide*.

- Caution: Sodium Azide is toxic.
 - All reagents contain Sodium Azide < 0,1% (w/v). Sodium Azide is not dangerous at this concentration. Careful handling is recommended to avoid accidental swallowing.
- Waste: Sodium Azide may form explosive compounds in metal drainage.

⚠️ Warnings

- All reagents are intended for *in vitro research*.
- Reagents have passed Quality Control checks and are expected to react as described above.
 - Shroud controls samples yield unexpected results, reagents are not to be considered reliable.
- Do not mix reagents from different lot numbers.
- Calibrators dilutions are stable for a few hours.

Technical Information and References

Please call Customer Support Service or visit www.newscientific.com.

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Technical Parameters

Calibration Curve Dilutions

Prepare various calibrator levels as shown in the table below:

Calibrator Dilutions (neat %)		2,5 %	5 %	10 %	20 %	40 %	Alternatively, prepare double the amount of the 40% dilution and attain 20%, 10%, 5%, and 2.5% doubling.
<input type="checkbox"/> Calibrator	µl	20	20	40	80	160	
<input type="checkbox"/> Buffered Saline (PBS)	µl	780	380	360	320	240	
Equivalent Concentrations	mg/dl	0,5	1	2	4	8	

To get set point 0 mg/dl, use saline.

Assay Parameters set up on APS

Program the assay parameters according to the following pattern

UPC ID	FRK (Free kappa) - FRL (Free lambda) or FKL (Free kappa+lambda - mixed)
Reagent position	Select free position
Gain	33 (for all reactions)
Optics	L (left) for FLC kappa - R (right) for FLC lambda, or vice versa L (left) or R (right) for FLC kappa+lambda - mixed
Dilution	Neat


Results Sample Rate - (minus) Reagent Blank rate

Programme the samples selecting the "User Prepared" Dilution.

It is a good idea to arrange a "Panel", dilution: "User Prepared"

Calibration Curve Examples

FRK				Date	FRL				Date	
				02/02/2001					02/02/2001	
Reagent	R.APS.FRK	Lot	1078-A012		Reagent	R.APS.FRL	Lot	1011-A102		
Calibrator	W.UPC.FRK	Lot	A012		Calibrator	W.UPC.FRL	Lot	A012		
Procedure		Standard	HS		Procedure		Standard	HS		
	mg/dl	rate	CV%			mg/dl	rate	CV%		
Saline	0	22			Saline	0	14			
Calibrator	0,25				Calibrator	0,25				
	0,5	64				0,5	124			
	1	145				1	362			
	2	354				2	834			
	4	717				4	1560			
	8	1530				8	2530			


Please Note: the value of the curves shown are only indicative

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New Scientific Company S.r.l. – Via Dante Alighieri, 35 – I 20032 Cormano (MI) – Italy