Italiano

Scopo

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English

Intended use

The Free Light Chains BNA Kit, when used in conjunction with Dade Behring "BN™ series" Analyzers (BNA, BN 100 BNII, ProSpec), is intended for the determination of *Free Light Chains (FLC) kappa and lambda (separately) in unconcentrated urine*.

Free Light Chains in Blood

By using the kits on serum and plasma, with the same application as that of urine, technically and analytically excellent results are obtained; however, we must emphasise that the clinical significance of the determination of the concentration of Free Light Chains in blood is currently being studied and defined. Thus this brochure refers only to urine.

For more information please call Customer Support Service or visit www.newscientific.com.

Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with absorbed antibodies specific to "hidden" sites of Light Chains of Immunoglobulins.

Contents

Kit Type					K
<i>Antiserum Reagent</i>					
<input type="checkbox"/> Free Light Chains kappa	<i>R.BNA.FRK</i>	green	1 vial – 4 ml	•	
<input type="checkbox"/> Free Light Chains lambda	<i>R.BNA.FRL</i>	yellow	1 vial – 4 ml	•	
<i>Calibrator</i>					
<input type="checkbox"/> Free Light Chains kappa	<i>W.UPC.FRK</i>	green	1 vial – 1 ml	•	
<input type="checkbox"/> Free Light Chains lambda	<i>W.UPC.FRL</i>	yellow	1 vial – 1 ml	•	
<i>Controls</i>					
<input type="checkbox"/> Free Light Chains kappa Control	<i>Y.UPR.FRK</i>	green	1 vial – 1 ml	-	
<input type="checkbox"/> Free Light Chains lambda Control	<i>Y.UPR.FRL</i>	yellow	1 vial – 1 ml	-	
• Reagents included in the kit.					
- Correlated reagents not included in the kit. They are available separately.					

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

Information re Calibrators / Controls

Due to lack of provisions from the Scientific Committees our "Internal Standardisation" has been applied.

See the relevant leaflet enclosed.

Analysr Set up

- BNA and BN100*** - Follow the manufacturer's instructions to program a new assay. The parameters are on the back.

- BNII and ProSpec*** - The protocol, which is analogous to that for "BNA", is for "pre-reaction" and, on BNII, for "interpolation at zero". For the parameters, contact Customer Services at NSC. For implementation, contact Dade Behring.

Processing Procedure

Calibration Procedure

The calibration procedure must be run for each chemistry, *The calibration curve must be rerun whenever a new lot number is used*. Results and curves on "BNII" are reported on the back.

Sample

Unconcentrated urine, centrifuged (3,000 rpm x 10'). Use 0,1% Sodium Azide as preservative.

Assay on samples

Samples must be run with FLC kappa and FLC lambda reagents.

Quality Control

For each set of determinations it is recommended to include:

<input type="checkbox"/> Negative Control	PBS (or Saline) = Reagent Blank.
<input type="checkbox"/> Positive Control	<input type="checkbox"/> Calibrator dilution 1:16 - Conc. ≈ 1.25 mg/dl. <p>The calibrator dilution is stable for a few hours.</p> <input type="checkbox"/> As a Positive Control alternative use: Y.UPR.FRK and Y.UPR.FRL

Results Interpretation

Reference Range

Free Light Chains are virtually absent in the urine of a healthy subject.

Significant Values.

Results higher than PBS results (Reagent Blank, Negative Control) are to be considered potentially positive.

Double PBS amount can be considered as cut-off. 1 mg/dl is generally assumed as cut-off for screening purposes.

⚠ ***Limitations – Lack of sample/calibrator parallelism***

The Free Light Chains in the sample could be "different" from those in calibrators and therefore they could react unproportionately; *lack of sample/calibrator parallelism*. The same occurs with the quantification of serum monoclonal immunoglobulins (serum Monoclonal Component – MC) by liquid phase ImmunoPrecipitation (turbidimetry, nephelometry).

⚠ ***Bit Data – Delta Bit***

Regular checking of the "Bit Data" - "delta bit" is recommended.

On BNA, BN100 and ProSpec, concentration below the calibration curve minimum (result < x mg/dl) can correspond to a "delta bit" just below the first point on the calibration curve; such a result could be significant.

⚠ ***Sample dilution***

Regular checking of the "Sample Dilution" is recommended.

⚠ ***Antigen Excess***

The application on "BNA" and "BN 100" allows a precise assay of a calibrator concentration ten times higher than the highest set point of the calibration curve. The application on "BNII" and "ProSpec" checks the antigen excess. However, Antigen Excess cannot be ruled out.

⚠ ***Monoclonality***

This test does not allow highlighting if the Free Light Chains in the sample are monoclonal or polyclonal.

⚠ Advice on Reagent usage

100 tests can be run on one reagent vial. Reagent vial dead volume is 100 µl.

⚠ Precautions and Hazards

Preservative: *Sodium Azide*.

- Caution: Sodium Azide is toxic.

All reagents contain Sodium Azide < 0,1% (w/v). *Sodium Azide is not dangerous at this concentration*. Careful handling is recommended to avoid accidental swallowing.
- Waste: Sodium Azide may form explosive compounds in metal drainage.

⚠ Warnings

- All reagents are intended for *in vitro research only*.
- Reagents have passed Quality Control checks and are expected to react as described above. Should controls samples yield unexpected results, reagents are not to be considered reliable.
- Do not mix reagents from different lot numbers.
- Calibrators dilutions are stable for a few hours.

Technical Information and References

Please call Customer Support Service or visit www.newscientific.com.

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Technical Parameters

Assay Parameters setup on BNA - BN100						
Reaction Volumes	Measurement Parameters			Calibration		
Sample	20	µl	Measure Type	Fixed Time	Set Points (Levels)	5
Buffer (*)	160	µl	First Measure	6 sec	First dilution	1:2,5
Reagent	40	µl	Last Measure	18 min	Sample dilution	1:1
(*) = N reaction buffer						

Assay Parameters setup on BNII - ProSpec							
Reaction Volumes	Measurement Parameters			Calibration			
Sample	pre reaction	reaction	Last Measure	pre reaction	reaction	Set Points (Levels)	
	6	µl	16	µl	180 sec	900 sec	5
Buffer (*)	200	µl	20	µl	Measure Type	Fixed Time	BNII only
Reagent	40	µl				extrapolation at "0" bit	
(*) = N reaction buffer							
Note	Complete BNII assay parameters are on document code:						
	NSC.FLC.AP.BNII	Valid revision	A31125				

FLC – Calibration Curve Example on BNII						
<i>Free Light Chains kappa</i>						
Analisi:133 Kf6 / curva di calibrazione di:25-11-2003 15:25:18 (curva attiva)						
Valore [Bit]						
Conc[mg/dL]						
Calibratore/Reagente	No.	Diluzione	Conc[mg/dL]	Val. preli. [Bit]	Valore iniziale[bit]	Valore [Bit]
68510 Kf6c standard	1		0,0000	0	0	0,00
	2	1: 40,0	0,5000	8	32198	-0,08
	3	1: 20,0	1,0000	18	445291	-1,78
6665 1 0 Free Kappa	4	1: 10,0	2,0000	25	362933	2,07
	5	1: 5,0	4,0000	93	4112281	-1,92
	6	1: 2,5	8,0000	259	5354463	0,63
<div style="display: flex; justify-content: space-between;"> BNII 2.2i Mar. 25 nov 2003, 15:43 </div>						
<i>Free Light Chains lambda</i>						
Analisi:154 lam6 / curva di calibrazione di:25-11-2003 15:33:02 (curva attiva)						
Valore [Bit]						
Conc[mg/dL]						
Calibratore/Reagente	No.	Diluzione	Conc[mg/dL]	Val. preli. [Bit]	Valore iniziale[bit]	Valore [Bit]
686010 Lf6c Standard	1		0,0000	0	0	0,00
	2	1: 40,0	0,5000	-1	186181	-0,28
	3	1: 20,0	1,0000	0	170328	0,77
6667 1 0 Lambda Free	4	1: 10,0	2,0000	13	1811783	-0,93
	5	1: 5,0	4,0000	165	2954064	0,76
	6	1: 2,5	8,0000	568	6705278	0,27
<div style="display: flex; justify-content: space-between;"> BNII 2.2i Mar. 25 nov 2003, 15:44 </div>						
Please note: the value of the curves are only indicative						

REF	K.BNA.FRK.FRL	IVD	
<p>Reagenti per uso professionale Reagents for professional use Reagenzien für Berufsgebrauch Reactivos para uso profesional Réactifs pour emploi professionnel Reagentes para uso profissional</p>			
<p>Kit Catene Leggere Libere – Reagenti separati Kit Free Light Chains – Separate Reagents Kit Freie Leichtketten – Separate Reagenzien Kit Cadenas Ligeras Libres – Reactivos Separados Kit Chaînes Légères Livres – Réactifs séparés Kit Cadeias Leves Livres – Reagentes Separados</p>			
<p>Istruzioni Instrucciones</p>			
<p>Instructions Instruccions</p>			
<p>Gebrauchsanweisung Instruções</p>			
<p>2 ... 8 °C</p>			
<p>100 Tests</p>			
<p>IT – EN – DE ES – FR – PT</p>			
<p>www.newscientific.com</p>			

REF	K.BNA.FRK.FRL	IVD	
<p>Reagenti per uso professionale Reagents for professional use Reagenzien für Berufsgebrauch Reactivos para uso profesional Réactifs pour emploi professionnel Reagentes para uso profissional</p>			

REF	K.BNA.FRK.FRL	IVD	
<p>Reagenti per uso professionale Reagents for professional use Reagenzien für Berufsgebrauch Reactivos para uso profesional Réactifs pour emploi professionnel Reagentes para uso profissional</p>			
<p>Kit Catene Leggere Libere – Reagenti separati Kit Free Light Chains – Separate Reagents Kit Freie Leichtketten – Separate Reagenzien Kit Cadenas Ligeras Libres – Reactivos Separados Kit Chaînes Légères Livres – Réactifs séparés Kit Cadeias Leves Livres – Reagentes Separados</p>			
<p>Istruzioni Instrucciones</p>			
<p>Instructions Instruccions</p>			
<p>Gebrauchsanweisung Instruções</p>			
<p>2 ... 8 °C</p>			
<p>100 Tests</p>			
<p>IT – EN – DE ES – FR – PT</p>			
<p>www.newscientific.com</p>			

REF	K.BNA.FRK.FRL	
------------	---------------	---------