Italiano

Scopo

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English

Indented Use

The Properdin factor B (PFB, GBG, C3PA) "BNA" kit , when used in conjunction with Dade Behring "BN™ series" Analyzers (BNA, BN100, BNII, ProSpec), is intended for the determination of Properdin factor B in serum.

Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with antibodies specific to Properdin factor B.

Reagents

Kit Type			KR	K
Antiserum Reagent				
<input type="checkbox"/> Properdin factor B	<i>R.BNA.GBG</i>	1 vial – 2 ml	•	•
Calibrator				
<input type="checkbox"/> Complement system protein	<i>W.SPC.CS1</i>	1 vial – 1 ml	-	•
Controls				
<input type="checkbox"/> Complement system proteins Medium Value	<i>Y.SPR.CS1.MV</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Complement system proteins Low Value	<i>Y.SPR.CS1.LV</i>	1 vial – 1 ml	-	•
• Reagents included in the kit.				
- Correlated reagents not included in the kit. They are available separately.				

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

Information re Calibrators / Controls

Due to lack of provisions from the Scientific Committees our "Internal Standardisation" has been applied.

See the relevant leaflet enclosed.

Analysar Set up

BNA and BN100 - Follow the manufacturer's instructions to program one new assay. The parameters are on the back.

BNII and ProSpec - The protocol is analogous to that for "BNA". For the parameters, contact Customer Services at NSC. For implementation, contact Dade Behring.

Processing Procedure

Calibration Curve

The calibration procedure must be run for each chemistry.

The calibration curve must be rerun whenever a new lot number is used.

Results and curves are reported on the back.

Sample

Serum: Fresh serum, limpid, not hemolyzed, not lipemic.

Assay on samples

Samples must be run as provided for by analyser instruction.

Quality Control

For each set of determinations it is recommended to include:

<input type="checkbox"/> Control Medium value	Y.SPR.CS1.MV
<input type="checkbox"/> Control Low value	Y.SPR.CS1.LV

Results Interpretation

Reference Range

It is advisable that each laboratory settles its own normal reference range.

The literature report: 10 – 45 mg/dl (Putnam F.W. – The Plasma Protein – II Edition).

⚠ Advice on Reagent usage

50 tests can be run on one reagent vial. Reagent vial dead volume is 100 µl.

⚠ Precautions and Hazards

Preservative: SodioAzide.

▪ Caution: Sodium Azide is toxic.

All reagents contain Sodium Azide < 0,1% (w/v). **Sodium Azide is not dangerous at this concentration.**

Careful handling is recommended to avoid accidental swallowing.

▪ Waste: Sodium Azide may form explosive compounds in metal drainage.

⚠ Warnings

▪ All reagents are intended for *in vitro* "*research*".

▪ Reagents have passed Quality Control checks and are expected to react as described above.

▪ Should controls samples yield unexpected results, reagents are not to be considered reliable.

▪ Do not mix reagents from different lot numbers.

▪ Calibrators dilutions are stable for a few hours.

Technical Information

Please call Customer Support Service or visit www.newscientific.com.

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New Scientific Company S.r.l. - Via Dante Alighieri, 35 - I 20032 Cormano (MI) - Italy