ESPAÑOL	(rev.: 09.01.2004)
---------	--------------------

OBJETO

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INTENDED USE

The *Anti-Streptolysin O - Kit DIME* Kit, when used with the *Dimension®* series Clinical Chemistry analyzers, of *Dade Behring*, is intended for the determinarion of **Anti-Streptolysin O (ASL)** in human serum.

METHOD

The method is based on the immunoprecipitation reaction in liquid phase (turbidimetry) enhanced with latex particles.

CONTENTS - KIT COMPONENTS

Latex Antiserum Reagent:

- Anti-Streptolysin O - Latex* (R.DIME.ASLAS) 2 vials of 3 ml

Reaction Buffer:

- Anti-Streptolysin O - Buffer* (R.DIME.ASL.BF) 2 vials of 18 ml

Form: líquido, ready for use.

Stability: 12 months at +2°C ... +8°C.

OTHER RELATED REAGENTS, NOT INCLUDED IN THE KIT

Calibrators:

- Anti-Streptolysin O - Calibrator* (W.ASL.STD) 3 vials of 1 ml

Controls:

- Anti-Streptolysin O - High Control* (W.ASL.CTH) 6 vials of 1 ml
- Anti-Streptolysin O - Low Control* (W.ASL.CTL) 6 vials of 1 ml

Information about Calibrators and Controls:

See the relevant Instructions.

ANALYZER SET UP

Follow the manufacturer of the analyzer instructions to set up a method (on Open Channels). Parameters can be found on the back

OPERATING PROCEDURE

Reagent Preparation - Filling the *Flex®*

Reagents supplied are enough to fill two *Flex®* (Reagent Cartridge) of 50 tests.

To fill a *Flex®* (see back), drain the content of a vial of Latex Antiserum Reagent (Comp. B) (R1-C2) in the well #1 and fill wells from #2 to #6 with 3.5 ml of Reaction Buffer (Comp. A) (R1-C1) (distribute the content of a vial between the five wells).

It is recommended to use a fine point pipette and proceed carefully to avoid reagent loss filling cartridges.

Calibration Curve

Set up the calibrator (W.ASL.STD) dilutions as shown in the table on the back.

Calibration Procedure

Run the calibration procedure preveded for the analyzer.

For the calibration curve construction use the "Logit" type (Not Lineal - 5 points).

The calibration must be rerun, at least, with each new lot.

The results and a tipycal curve are reported as an example on the back.

Samples

Fresh Serum, limpid, not hemolyzed, not lipemic.

Assay on Samples

Samples must be run as preveded for the analyzer.

Quality Control

It is advisable to insert in each analytical serie:

- High Value Control: W.ASL.CTH
- Low Value Control: W.ASL.CTL

RESULTS AND INTERPRETATION

Reference Values

It is advisable that each laboratory settles its own normal reference range.

ASL reference values specified in literature are up to 200 IU/ml (300 IU/ml for children).

Antigen Excess

The proposed application allows an accurate assay of a calibrator concentration four times higher than the highest set point of the calibration curve (2800 IU/ml). However, "Antigen Excess" cannot be ruled out.

ADVICE ON REAGENT CONSUMPTION

The volume of reagents supplied in the vials allows to run 100 test with the proposed standard parameters; to avoid loss filling the reagent containers it is recommended to use a fine point pipette.

PRECAUTIONS AND HAZARDS

- Preservative:** *Sodium Azide*.
 - Handling:** Sodium Azide is toxic. All the reagents contain Sodium Azide <0,1% (w/v). Even if Sodium Azide is not dangerous at this concentration, make sure to take the necessary precautions to prevent accidental swallowing.
 - Waste:** Sodium Azide can react with lead and copper plumbing to form highly explosive metal azides.
- Since absence of infectious agents can not be proven with absolute certainty, all materials derived from human or animal blood should always be handled with due care, following the precautions recommended for biohazardous material.

WARNINGS

- All reagents are intended only for "*In Vitro* research".
- Reagents have passed quality control cheks and are expected to react as described in this Instructions. If controls do not yield the preveded reaction, all the reagents would not be considered reliable.
- Do not mix reagents from different lot kits.
- Calibrators dilutions are stable for a few hours.

TECHNICAL INFORMATION

Contact Customer Support Service or visit www.newscientific.com.

PORTUGUÊS	(rev.: 23.06.2004)
-----------	--------------------

PORTUGUÊS	(rev.: 23.06.2004)
-----------	--------------------

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CARTRIDGE CONFIGURATION (FLEX® FILLING)

well #1 Component B 3.0 ml Anti-Streptolysin O - Latex (REF: R.DIME.ASL.AS)
wells #2 to #6 Component A 3.5 ml Anti-Streptolysin O - Buffer (REF: R.DIME.ASL.BF)

ASSAY PARAMETERS SET UP ON DIMENSION®

User-Defined Method


CHANNEL: User selected NAME: User defined (x_...)
MEASUREMENT MODE: **TURBIDIMETRIC** STANDARD CURVE: **LOGIT**
SAMPLE (S1): TIME: **0.0** sec CHASE: **16** ul Mix: **STRONG**
VOLUME: **4** ul
FIRST REAGENT (R1): TIME: **-60.0** sec CHASE: **26** ul Mix: **STRONG**
COMPONENT 1: (A) **300** ul COMPONENT 2: (B) **54** ul COMPONENT 3: () **0** ul
SECOND REAGENT (R2): not defined
THIRD REAGENT (R3): not defined
PHOTOMETRY TIMES: P1: **-10.0** sec P2: **400.0** sec P3: not defined P4: not defined
REAGENT CARTRIDGE: WELL: 1 2 3 4 5 6
COMPONENT: (B) (A) (A) (A) (A) (A)
ALIQUOTS: **50** **10** **10** **10** **10** **10**
LIFE (HRS): **720** **720** **720** **720** **720** **720**
REAG. CART. LIFE: **720** hrs CALIBR. INTERVAL: user defined (maximum 2160 hrs)
MAU CALCULATION: MODE: **FINAL POINT** MAIN FILTER: **540** SECOND. FILTER: not defined
P1 (-10 SEC): DIL. FACTOR: **-0.950** DOI: **0.00 (INACTIVE)**
P2 (400 SEC): DIL. FACTOR: **1.000** DOF: **0.00 (INACTIVE)**

Method Parameters

NAME: User defined DECIMALS: User defined UNITS: User defined CALCULATION: **LOGIT**
VOL. SELEC.: **4** ul VOL. AUTODIL.: **2** ul VOL. URINE: **0** ul
INTERVALS: REFERENCE: User defined (Warning! forced field definition)
ASSAY: User defined (recommended from 90 to 700 IU/ml (#1 and #5 calibration point values)) (Warning! forced field definition)
PANIC: User defined (usually the same as Assay Range)
COEF. (C0 TO C4): Calculated by the analyzer (for the first time method calibration set coefficients to C0=1, C1=2, C2=3, C3=4 and C4=0.5)

EXAMPLE OF CALIBRATION CURVE

REAGENTS: REF.: K.DIME.ASL - LOT: 3362 CALIBRATOR: REF.: W.ASL.STD - LOT: 3383

Note - Warning

The results of the calibration and the graphic have value only as an example

BIOMEDICINA - PATOLOGIA CLINICA - INFORMATICA

Immuno Turbidimetric System - on **DIMENSION®**

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New Scientific Company España, S.r.l.