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### ASSAY PARAMETERS SET UP ON DIMENSION®

#### User-Defined Method


CHANNEL: User selected NAME: User defined (x...)  
MEASUREMENT MODE: TURBIDIMETRIC STANDARD CURVE: LOGIT  
SAMPLE (S1): TIME: 0.0 sec CHASE: 10 ul Mix: STRONG  
VOLUME: 20 ul  
FIRST REAGENT (R1): TIME: -60.0 sec CHASE: 0 ul Mix: NONE  
COMPONENT 1: (A) 320 ul COMPONENT 2: ( ) 0 ul COMPONENT 3: ( ) 0 ul  
SECOND REAGENT (R2): TIME: 90.0 sec CHASE: 30 ul Mix: STRONG  
COMPONENT 1: (B) 70 ul COMPONENT 2: ( ) 0 ul COMPONENT 3: ( ) 0 ul  
PHOTOMETRY TIMES: P1: 80.0 sec P2: 300.0 sec P3: not defined P4: not defined  
REAGENT CARTRIDGE: WELL: 1 2 3 4 5 6  
COMPONENT: (B) (A) (A) (A) (A) (A)  
ALIQUOTS: 50 10 10 10 10 10  
LIFE (HRS): 720 720 720 720 720 720  
REAG. CART. LIFE: 720 hrs CALIBR. INTERVAL: user defined (maximum 2160 hrs)  
MAU CALCULATION: MODE: FINAL POINT MAIN FILTER: 452 SECOND. FILTER: 700  
P1 (80 SEC): DIL. FACTOR: -0.777 DOI: 0.00 (INACTIVE)  
P2 (300 SEC): DIL. FACTOR: 1.000 DOF: 0.00 (INACTIVE)

#### Method Parameters

NAME: User defined DECIMALS: User defined UNITS: User defined CALCULATION: LOGIT  
VOL. SELEC.: 20 ul VOL. AUTODIL.: 5 ul VOL. URINE: 5 ul  
INTERVALS: REFERENCE: User defined (Warning! forced field definition)  
ASSAY: User defined (#1 and #5 calibration point values (the values are calibrator lot dependent)) (Warning! forced field definition)  
PANIC: User defined (usually the same as Assay Range)  
COEF. (C0 to C4): Calculated by the analyzer (for the first time method calibration set coefficients to C0=1, C1=2, C2=3, C3=4 and C4=0.5)

### EXAMPLE OF CALIBRATION CURVE

REAGENTS: REF.: K.DIME.MALB - LOT: 3331 CALIBRATOR: REF.: W.MALB.STD - LOT: 3314

### Note - Warning

The results of the calibration and the graphic have value only as an example

BIOMEDICINA - PATOLOGIA CLINICA - INFORMATICA

Immuno Turbidimetric System - on DIMENSION®

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**INTENDED USE**

The *Micro Albumin (urine) - Kit DIME* Kit, when used with the *Dimension®* series Clinical Chemistry analyzers, of *Dade Behring*, is intended for the determination of **Micro Albumin (mALB)** in human urine.

**METHOD**

The method is based on the immunoprecipitation reaction in liquid phase (turbidimetry) with antibodies specific to the Albumin.

**CONTENTS - KIT COMPONENTS**

- Antiserum Reagent:**
- *Micro Albumin - Antiserum* (R.DIME.MALB.AS) 2 vials of 4 ml
- Reaction Buffer:**
- *Micro Albumin - Buffer* (R.DIME.MALB.BF) 2 vials of 18 ml
- Form:** liquid, ready for use.  
**Stability:** 12 months at +2°C ... +8°C.

**OTHER RELATED REAGENTS, NOT INCLUDED IN THE KIT**

- Calibrators:**
- *Micro Albumin (urine) - Calibrator* (W.MALB.STD) 3 vials of 1 ml
- Controls:**
- *Micro Albumin (urine) - High Control* (W.MALB.CTH) 6 vials of 1 ml
  - *Micro Albumin (urine) - Low Control* (W.MALB.CTL) 6 vials of 1 ml

**Information about Calibrators and Controls:**

See the relevant Instructions.

**ANALYZER SET UP**

Follow the manufacturer of the analyzer instructions to set up a method (on Open Channels). Parameters can be found on the back

**OPERATING PROCEDURE**

**Reagent Preparation - Filling the *Flex®***

Reagents supplied are enough to fill two *Flex®* (Reagent Cartridge) of 50 tests. To fill a *Flex®* (see back), drain the content of a vial of Antiserum Reagent (Comp. B) (R2-C1) in the well #1 and fill wells from #2 to #6 with 3.5 ml of Reaction Buffer (Comp. A) (R1-C1) (distribute the content of a vial between the five wells). It is recommended to use a fine point pipette and proceed carefully to avoid reagent loss filling cartridges.

**Calibration Curve**

Set up the calibrator (W.MALB.STD) dilutions as shown in the table on the back.

**Calibration Procedure**

Run the calibration procedure previded for the analyzer. For the calibration curve construction use the "Logit" type (Not Lineal - 5 points). **The calibration must be rerun, at least, with each new lot.**

The results and a typical curve are reported as an example on the back.

**Samples**

Fresh Urine, limpid, centrifuged (3000 rpm - 10')..

**Assay on Samples**

Samples must be run as previded for the analyzer.

**Quality Control**

It is advisable to insert in each analytical serie:

- High Value Control: W.MALB.CTH
- Low Value Control: W.MALB.CTL

**RESULTS AND INTERPRETATION**

**Reference Values**

It is advisable that each laboratory settles its own normal reference range. mALB reference values specified in literature are variable, depending on the authors and the clinical criteria followed,

and usually are expressed in relation to the diuresis (mg/24h, µg/min or referred to the Creatinine excretion). An excretion higher than 20 mg/24h is usually considered pathologic.

**Antigen Excess**

The proposed application allows an accurate assay of a calibrator concentration four times higher than the highest set point of the calibration curve (800 - 1000 mg/l). However, "Antigen Excess" cannot be ruled out.

**ADVICE ON REAGENT CONSUMPTION**

The volume of reagents supplied in the vials allows to run 100 test with the proposed standard parameters; to avoid loss filling the reagent containers it is recommended to use a fine point pipette.

**PRECAUTIONS AND HAZARDS**

- **Preservative:** *Sodium Azide*.
  - **Handling:** Sodium Azide is toxic. All the reagents contain Sodium Azide <0,1% (w/v). Even if Sodium Azide is not dangerous at this concentration, make sure to take the necessary precautions to prevent accidental swallowing.
  - **Waste:** Sodium Azide can react with lead and copper plumbing to form highly explosive metal azides.
- Since absence of infectious agents can not be proven with absolute certainty, all materials derived from human or animal blood should always be handled with due care, following the precautions recommended for biohazardous material.

**WARNINGS**

- All reagents are intended only for "In Vitro research".
- Reagents have passed quality control cheks and are expected to react as described in this Instructions. If controls do not yield the previded reaction, all the reagents would not be considered reliable.
- Do not mix reagents from different lot kits.
- Calibrators dilutions are stable for a few hours.

**TECHNICAL INFORMATION**

Contact Customer Support Service or visit www.newscientific.com.

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