Italiano

Scopo

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English

Indented Use

The C1q (C1Q) *IMG* Kit, when used in conjunction with Beckman Coulter Immage® Analyzer, is intended for the determination of C1q in serum.

Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with antibodies specific to C1q.

Kit Type		KR	K	
<i>Antiserum Reagent</i>				
<input type="checkbox"/> C1q	<i>R.IMG.C1Q</i>	1 vial – 3,5 ml	•	•
<i>Calibrator</i>				
<input type="checkbox"/> Complement system protein	<i>W.SPC.CS1</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Prediluted Calibrator Set <p>Complement system proteins A Set of ready-to-use Prediluted Calibrators is available.</p>	<i>W.SET.SPC.CS1</i>	5 vials – 1 ml	-	-
<i>Controls</i>				
<input type="checkbox"/> Complement system proteins <p>Medium Value</p>	<i>Y.SPR.CS1.MV</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Complement system proteins <p>Low Value</p>	<i>Y.SPR.CS1.LV</i>	1 vial – 1 ml	-	•
• Reagents included in the kit.				
- Correlated reagents not included in the kit. They are available separately.				

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

Information re Calibrators / Controls

Due to lack of previsions from the Scientific Committees our “Internal Standardisation” has been applied.

See the relevant leaflet enclosed.

Analysér Set up

Follow the manufacturer’s instructions to program one “UDR” assay. The assay parameters can be found on the back.

Getting the cartridges ready

One UDR cartridge must be labelled and filled as above:

- Compartment A: with Antiserum Reagents (*R.IMG.C1Q*);

Processing Procedure

Calibration Curve

Prepare various calibrator levels (see table): calibrators dilutions are stable for a few hours.

As an alternative use the Prediluted Calibrator Set – W.SET.SPC.CS1.

Calibration Procedure

The calibration procedure must be run for the chemistry.

The calibration curve must be rerun at every cartridge replacement and whenever a new lot number is used.

Results and curves are reported on the back.

Sample

- Serum: Fresh serum, limpид, not hemolyzed, not lipemic.

Assay on samples

Samples must be run as provided for by analyser instruction.

Quality Control

For each set of determinations it is recommended to include:

<input type="checkbox"/> Control Medium value	Y.SPR.CS1.MV
<input type="checkbox"/> Control Low value	Y.SPR.CS1.LV

Results Interpretation

Reference Range

It is advisable that each laboratory settles its own normal reference range.

The literature report: 10 – 25 mg/dl (Putnam F.W. – The Plasma Protein – II Edition).

⚠️ Advice on Reagent usage

100 tests can be run on one reagent vial; please use serological pipette to avoid reagent leaks filling cartridges.

Cartridges Compartment dead volume are 250 µl.

⚠️ Precautions and Hazards

Preservative: SodioAzide.

- Caution: Sodium Azide is toxic.
 - All reagents contain Sodium Azide < 0,1% (w/v). **Sodium Azide is not dangerous at this concentration.** Careful handling is recommended to avoid accidental swallowing.
 - Waste: Sodium Azide may form explosive compounds in metal drainage.

⚠️ Warnings

- All reagents are intended for ***in vitro "research"***.
- Reagents have passed Quality Control checks and are expected to react as described above. Should controls samples yield unexpected results, reagents are not to be considered reliable.
- Do not mix reagents from different lot numbers.
- Calibrators dilutions are stable for a few hours.

Technical Information

Please call Customer Support Service or visit www.newsscientific.com.

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Technical Parameters

C1Q - Calibration Curve Dilution						
REF	W.SPC.CS1	LOT	<i>A501</i>			

Prepare various calibrator levels doubling or as shown in the table below:

Calibrator Dilutions	1 : 16	1 : 8	1 : 4	1 : 2	1 : 1	
<input type="checkbox"/> Calibrator	µl	25	50	50	100	200
<input type="checkbox"/> Buffered Saline (PBS)	µl	375	350	150	100	0

Equivalent Calibrator Concentration - Neat mg/dl (*)	36,0	2,3	4,5	9,0	18,0	36,0
--	------	-----	-----	-----	------	------

* The calibrator values are lot dependent.


As an alternative use Prediluted Calibrator Set: Code W.SET.SPC.CS1

To get set point 0 mg/dl use saline or Diluent 1

C1Q - Assay Parameters set up on Immage						
Protocol Information						
Lot Number	see "UDR" cartridge	Unit	mg/dl			
Reagent Serial	see "UDR" cartridge	Protocol	Non-Competitive Nephelometric			
Reagent Expiry Date	see Reagent vial					
Sample or Dilution Volume	50 µl	Gain	4			
Reaction Buffer Volume (1)	195 µl	Calibrator Dilution	1 : 5			
Compartment A Volume	30 µl	Sample Dilution	1 : 5 (2)			
Compartment B Volume	0 µl	Reaction Time	2,5 minutes			
		Notice:				
Levels	6	- Model	3rd order polynomial			
Replicates	2	- Set point 1 - Input	0,001			
(1) Reaction Buffer = Buffer 1		(2) To be configured after calibration approval				

C1Q - Calibration Curve Example

<i>C1Q</i>		Date		24/03/05		
Reagent	R.IMG.C1Q	Lot	A501			
Calibrator	W.SPC.CS1	Lot	A501			
Cal. Concentration	Curve	signal media	<i>CV%</i>	Sample mg/dl		
Net mg/dl	36	mg/dl				
	Saline	0,001	1	<i>0,28</i>	0,001	
Calibrator	1 : 16	2,25	132	<i>2,39</i>	2,25	
	1 : 8	4,50	341	<i>0,16</i>	4,50	
	1 : 4	9,00	839	<i>0,10</i>	9,00	
	1 : 2	18,00	1897	<i>0,23</i>	18,00	
	1 : 1	36,00	3818	<i>1,44</i>	36,00	






Please Note: the value of the curves shown are only indicative

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New Scientific Company S.r.l. – Via Dante Alighieri, 35 – I 20032 Cornoano (MI) – Italy