## Italiano

### Scopo

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## English

### Indented Use

The Free Light Chains "IMG" Kit, when used in conjunction with Beckman-Coulter Immage® Analyzer, is intended for the determination of *Free Light Chains (FLC) kappa and lambda (separately) in unconcentrated urine*.

*Free Light Chains in Blood*

By using the kits on serum and plasma, with the same application as that of urine, technically and analytically excellent results are obtained; however, we must emphasise that the clinical significance of the determination of the concentration of Free Light Chains in blood is currently being studied and defined. Thus this brochure refers only to urine.

For more information please call Customer Support Service or visit www.newscientific.com.

### Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with absorbed antibodies specific to "hidden" sites of Light Chains of Immunoglobulins.

### Contents

Kit Type				K
<i>Antiserum Reagents</i>				
<input type="checkbox"/> Free Light Chains kappa	<i>R.IMG.FRK</i>	green	1 vial – 2,4 ml	•
<input type="checkbox"/> Free Light Chains lambda	<i>R.IMG.FRL</i>	yellow	1 vial – 2,4 ml	•
<i>Accessory Reagent</i>				
<input type="checkbox"/> Reagent "B" FLC	<i>R.IMG.FLC.B</i>	azure	1 vial – 3 ml	•
<i>Calibrators</i>				
<input type="checkbox"/> Free Light Chains kappa	<i>W.UPC.FRK</i>	green	1 vial – 1 ml	•
<input type="checkbox"/> Free Light Chains lambda	<i>W.UPC.FRL</i>	yellow	1 vial – 1 ml	•
<input type="checkbox"/> Prediluted Calibrator Set				
<input type="checkbox"/> Free Light Chains kappa – lambda	<i>W.SET.UPC.FRK.FRL</i>		5 vials – 1 ml	-
<input type="checkbox"/> Set of ready-to-use Calibrators			5 vials – 1 ml	-
<i>Controls</i>				
<input type="checkbox"/> Free Light Chains kappa Control	<i>Y.UPR.FRK</i>	green	1 vial – 1 ml	-
<input type="checkbox"/> Free Light Chains lambda Control	<i>Y.UPR.FRL</i>	yellow	1 vial – 1 ml	-
<input type="checkbox"/> Reagents included in the kit.				
<input type="checkbox"/> - Correlated reagents not included in the kit. They are available separately.				

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

### Information re Calibrators / Controls

Due to lack of provisions from the Scientific Committees our "Internal Standardisation" has been applied.

See the relevant leaflet enclosed.

### Analysar Set up

Follow the manufacturer's instructions to program two "UDR" assays.

The assay parameters can be found on the back.

*Getting the cartridges ready*

Two UDR cartridges must be labelled and filled as above:

- Compartment A: with Antiserum Reagents;
- Compartment B: with about 1 ml of Accessory Reagent "B".

### Processing Procedure

Calibration Curve

Prepare various calibrator levels (see table) or use the Prediluted Calibrator Sets: W.SET.UPC.FRK.FRL.

*Calibration Procedure*

The calibration procedure must be run for each chemistry. Results and curves are reported on the back.

*The calibration curve must be rerun at every cartridge replacement and whenever a new lot number is used.*

*Sample*

Unconcentrated urine, centrifuged (3,000 rpm x 10'). Use 0,1% Sodium Azide as preservative.

### Assay on samples

Samples must be run with FLC kappa and FLC lambda reagents.

*Quality Control*

For each set of determinations it is recommended to include:

<input type="checkbox"/> Negative Control	PBS (or Saline) = Reagent Blank. <ul style="list-style-type: none"><li><input type="checkbox"/> Calibrator dilution 1:16 – Conc. ≈ 1.25 mg/dl.</li></ul>
<input type="checkbox"/> Positive Control	The calibrator dilution is stable for a few hours. <ul style="list-style-type: none"><li><input type="checkbox"/> As a Positive Control alternative use: Y.UPR.FRK and Y.UPR.FRL.</li></ul>

### Results Interpretation

*Reference Range*

Free Light Chains are virtually absent in the urine of a healthy subject.

*Significant Values*

Results higher than PBS results (reagent blank, negative control) are to be considered potentially positive. Double PBS amount can be considered as cut-off. 1 mg/dl is generally assumed as cut-off for screening purposes.

*⚠ Limitations – lack of sample/calibrator parallelism*

The Free Light Chains in the sample could be "different" from those in calibrators and therefore they could react unproportionately; *lack of sample/calibrator parallelism*. The same occurs with the quantification of serum monoclonal immunoglobulins (serum Monoclonal Component – MC) by liquid phase ImmunoPrecipitation (turbidimetry, nephelometry).

*⚠ Antigen Excess*

The application allows a precise assay of a calibrator concentration ten times higher than the highest set point of the calibration curve. However, Antigen Excess cannot be ruled out.

*⚠ Monoclonality*

This test does not allow highlighting if the Free Light Chains in the sample are monoclonal or polyclonal.

### ⚠ Advice on Reagent usage

100 tests can be run on one reagent vial; please use serological pipette to avoid reagent leaks filling cartridges.

Cartridges Compartment dead volume are 250 µl.

### ⚠ Precautions and Hazards

Preservative: *Sodium Azide*.

- Caution: Sodium Azide is toxic.
  - All reagents contain Sodium Azide < 0,1% (w/v). **Sodium Azide is not dangerous at this concentration**. Careful handling is recommended to avoid accidental swallowing.
- Waste: Sodium Azide may form explosive compounds in metal drainage.

### ⚠ Warnings

- All reagents are intended for *in vitro research*.
- Reagents have passed Quality Control checks and are expected to react as described above.
- Should controls samples yield unexpected results, reagents are not to be considered reliable.
- Do not mix reagents from different lot numbers.
- Calibrators dilutions are stable for a few hours.

### Technical Information and References

Please call Customer Support Service or visit www.newscientific.com.

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## Technical Parameters

FLC - Calibration Curve Dilution						
<b>REF</b>	W.UPC.FRK - W.UPC.FRL - W.UPC.FKL	<b>LOT</b>	A411			
Prepare various calibrator levels as shown in the table below:						
Calibrator Dilutions (neat 100%)		2,5%	5%	10%	20%	40%
		1:40	1:20	1:10	1:5	1:2,5
<input type="checkbox"/> Calibrator	µl	10	20	20	50	100
<input type="checkbox"/> Buffered Saline (PBS)	µl	390	380	180	200	150
Equivalent Calibrator Concentration						
Neat mg/dl (*)	20,0	0,5	1,0	2,0	4,0	8,0
* The calibrator values are lot dependent.						
Alternatively, prepare double the amount of the 40% dilution and attain 20%, 10%, 5%, and 2.5% doubling.						
As an alternative use Prediluted Calibrator Set, Code: W.SET.UPC.FRK.FRL						
To get set point 0 mg/dl use saline or Diluent 1						

FLC - Assay Parameters set up on Immage
---

## Protocol Information


Lot Number	see "UDR" cartridge	Unit	mg/dl
Reagent Serial	see "UDR" cartridge	Protocol	Non-Competitive
Reagent Expiry Date	see Reagent vial		Nephelometric

Method	Standard		
Sample or Dilution Volume	12 µl	Gain	3
Reaction Buffer Volume	195 µl	Calibrator Dilution	1:1
Compartment A Volume	20 µl	Sample Dilution	1:1
Compartment B Volume	5 µl	Reaction Time	2,5 minutes


## Calibration Information

Levels	6	Notice - Concentration Input	
Replicates	2	Lowest level, Set Point 1 - Input	0.0010

FLC - Calibration Curve Examples						
<i>FLC kappa</i>		Date	24/03/05			
Reagent	R.IMG.FRK	Lot	IMFK-028A501			
Calibrator	W.SET.FRK	Lot	A501			
Procedure	Standard	HS				
Saline	mg/dl	signal	CV%	signal	CV%	
	0	246	3,15			
Calibrator	0,25	296	1,40			
	0,5	341	1,16			
	1	455	0,58			
	2	655	0,72			
	4	1074	11,92			
	8	2130	2,03			
<i>FLC lambda</i>		Date	24/03/05			
Reagent	R.IMG.FRL	Lot	IMKL-A501			
Calibrator	W.SET.FRL	Lot	A501			
Procedure	Standard	HS				
Saline	mg/dl	signal	CV%	signal	CV%	
	0	111	10,67			
Calibrator	0,25	164	2,30			
	0,5	209	0,57			
	1	314	1,29			
	2	513	0,11			
	4	920	4,91			
	8	1539	0,10			

Please Note: the value of the curves shown are only indicative

BIOMICEDINA – PATOLOGIA CLINICA - INFORMATICA

<b>REF</b>	K.IMG.FRK.FRL	<b>IVD</b>	
			Reagenti per uso professionale
			Reagents for professional use
			Reagenzien für Berufsgebrauch
			Reactivos para uso profesional
			Réactifs pour emploi professionnel
			Reagentes para uso profissional