Italiano

Scopo

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English

Indented Use

The Free Light Chains "IMG" Kit, when used in conjunction with Beckman-Coulter Immage® Analyzer, is intended for the determination of *Free Light Chains (FLC) kappa and lambda (separately) in unconcentrated urine*.

Free Light Chains in Blood and Cerebral Spinal Fluid

By using the kits on serum (or plasma) and Cerebral Spinal Fluid, with the specific application, technically and analytically excellent results are obtained; however, we must emphasise that the clinical significance of the determination of the concentration of Free Light Chains in blood and cerebral Spinal Fluid is currently being studied and defined. Thus this brochure refers only to urine.

For more information please call Customer Support Service or visit www.newscientific.com.

Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with absorbed antibodies specific to "hidden" sites of Light Chains of Immunoglobulins.

Contents

Kit Type					K
<i>Antiserum Reagents</i>					
<input type="checkbox"/> Free Light Chains kappa	<i>R.IMG.FRK</i>	green	1 vial – 2,4 ml	●	
<input type="checkbox"/> Free Light Chains lambda	<i>R.IMG.FRL</i>	yellow	1 vial – 2,4 ml	●	
<i>Accessory Reagent</i>					
<input type="checkbox"/> Reagent "B" FLC	<i>R.IMG.FLC.B</i>	azure	1 vial – 3 ml	●	
<i>Calibrators</i>					
<input type="checkbox"/> Free Light Chains kappa	<i>W.UPC.FRK</i>	green	1 vial – 1 ml	●	
<input type="checkbox"/> Free Light Chains lambda	<i>W.UPC.FRL</i>	yellow	1 vial – 1 ml	●	
<input type="checkbox"/> Prediluted Calibrator Set <p>Free Light Chains kappa – lambda Set of ready-to-use Calibrators</p>	<i>W.SET.UPC.FRK.FRL</i>		5 vials – 1 ml <p>5 vials – 1 ml</p>	—	
<i>Controls</i>					
<input type="checkbox"/> Free Light Chains kappa Control - medium	<i>Y.UPR.FRK.MV-3</i>	green	1 vial – 3 ml	—	
<input type="checkbox"/> Free Light Chains kappa Control - low	<i>Y.UPR.FRK.LV-3</i>	green	1 vial – 3 ml	—	
<input type="checkbox"/> Free Light Chains lambda Control - medium	<i>Y.UPR.FRL.MV-3</i>	yellow	1 vial – 3 ml	—	
<input type="checkbox"/> Free Light Chains lambda Control - low	<i>Y.UPR.FRL.LV-3</i>	yellow	1 vial – 3 ml	—	
<input type="checkbox"/> Free Light Chains kappa Control	<i>Y.UPR.FRK</i>	green	1 vial – 1 ml	—	
<input type="checkbox"/> Free Light Chains lambda Control	<i>Y.UPR.FRL</i>	yellow	1 vial – 1 ml	—	
<input type="checkbox"/> Reagents included in the kit.					
— Correlated reagents not included in the kit. They are available separately.					

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

Information re Calibrators / Controls

Due to lack of precautions from the Scientific Committees our "Internal Standardisation" has been applied.

See the relevant leaflet enclosed.

Analysér Set up

Follow the manufacturer's instructions to program two "UDR" assays.

The assay parameters can be found on the back.

Getting the cartridges ready

Two UDR cartridges must be labelled and filled as above:

- Compartment A: with Antiserum Reagents;
- Compartment B: with about 1 ml of Accessory Reagent "B".

Processing Procedure

Calibration Curve

Prepare various calibrator levels (see table) or use the Prediluted Calibrator Sets: W.SET.UPC.FRK.FRL.

Calibration Procedure

The calibration procedure must be run for each chemistry. Results and curves are reported on the back.

The calibration curve must be rerun at every cartridge replacement and whenever a new lot number is used.

Sample

Unconcentrated urine, centrifuged (3,000 rpm x 10'). Use 0,1% Sodium Azide as preservative.

Assay on samples

Samples must be run with FLC kappa and FLC lambda reagents.

Quality Control

For each set of determinations it is recommended to include:

<input type="checkbox"/> Negative Control	PBS (or Saline) = Reagent Blank
<input type="checkbox"/> Positive Control	<input type="checkbox"/> Calibrator dilution 1:6 – Conc. ≈ 1.33 mg/dl. <p>The calibrator dilution is stable for a few hours.</p> <input type="checkbox"/> As a Positive Control alternative use: Y.UPR.FRK, Y.UPR.FRL; Y.UPR.FRK.MV-3, Y.UPR.FRL.MV-3; Y.UPR.FRK.LV-3, Y.UPR.FRL.LV-3

Results Interpretation

Reference Range

Free Light Chains are virtually absent in the urine of a healthy subject.

Significant Values

Results higher than PBS results (reagent blank, negative control) are to be considered potentially positive. Double PBS amount can be considered as cut-off. 1 mg/dl is generally assumed as cut-off for screening purposes.

⚠ Limitations – lack of sample/calibrator parallelism

The Free Light Chains in the sample could be "different" from those in calibrators and therefore they could react unproportionately; *lack of sample/calibrator parallelism*. The same occurs with the quantification of serum monoclonal immunoglobulins (serum Monoclonal Component – MC) by liquid phase ImmunoPrecipitation (turbidimetry, nephelometry).

⚠ Antigen Excess

The application allows a precise assay of a calibrator concentration ten times higher than the highest set point of the calibration curve. However, Antigen Excess cannot be ruled out.

⚠ Monoclonality

This test does not allow highlighting if the Free Light Chains in the sample are monoclonal or polyclonal.

⚠ Advice on Reagent usage

100 tests can be run on one reagent vial; please use serological pipette to avoid reagent leaks filling cartridges.

Cartridges Compartment dead volume are 250 µl.

⚠ Precautions and Hazards

Preservative: *Sodium Azide*.

- Caution: Sodium Azide is toxic.
 - All reagents contain Sodium Azide < 0,1% (w/v). **Sodium Azide is not dangerous at this concentration**. Careful handling is recommended to avoid accidental swallowing.
- Waste: Sodium Azide may form explosive compounds in metal drainage.

⚠ Warnings

- All reagents are intended for *in vitro research*.
- Reagents have passed Quality Control checks and are expected to react as described above. Should controls samples yield unexpected results, reagents are not to be considered reliable.
- Do not mix reagents from different lot numbers.
- Calibrators dilutions are stable for a few hours.

Technical Information and References

Please call Customer Support Service or visit www.newscientific.com.

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Technical Parameters

FLC - Calibration Curve Dilution - Urine							
REF	W.UPC.FRK - W.UPC.FRL - W.UPC.FKL	LOT	WFKA603 - WFKA603 - WFKLA603				
Prepare various calibrator levels as shown in the table below:							
Calibrator Dilutions (neat 100%)		6,25%	12,5%	25%	50%	100%	
		1:16	1:8	1:4	1:2	1:1	
<input type="checkbox"/> Calibrator	µl	25	50	100	200	400	
<input type="checkbox"/> Buffered Saline (PBS)	µl	375	350	300	200	0	
Equivalent Calibrator Concentration Neat mg/dl (*)		8,0	0,5	1,0	2,0	4,0	8,0
* The calibrator values are lot dependent. Alternatively, prepare double the amount of the 50% dilution and attain 25%, 12.5%, 6.25% doubling. As an alternative use Prediluted Calibrator Set, Code: W.SET.UPC.FRK.FRL To get set point 0 mg/dl use saline or Diluent 1							

FLC - Assay Parameters set up on Image - Urine					
Protocol Information					
Lot Number	see "UDR" cartridge	Unit	mg/dl		
Reagent Serial	see "UDR" cartridge	Protocol	Non-Competitive Nephelometric		
Reagent Expiry Date	see Reagent vial				
Method	Standard				
Sample or Dilution Volume	12 µl	Gain	3		
Reaction Buffer Volume (1)	195 µl	Calibrator Dilution	1 : 1		
Compartment A Volume	20 µl	Sample Dilution	1 : 1		
Compartment B Volume	5 µl	Reaction Time	2,5 minutes		
Calibration Information					
Levels	6		Notice:		
Replicates	2		- Set point 1 - Input 0,001		
			- Model 3 rd order polynomial		
(1) Reaction Buffer = Buffer 1					

FLC - Calibration Curve Examples - Urine					
<i>FLC kappa</i>			Date	10/04/06	
Reagent	R.IMG.FRK	Lot	IMFK-A603		
Calibrator	W.UPC.FRK	Lot	WFK-A603		
	Procedure	Standard	HS		
	mg/dl	signal	CV%	signal	CV%
	Saline	0	86	3,62	
		0,125	107	4,01	
		0,250	120	3,25	
		0,50	148	2,46	
		1,00	216	0,52	
		2,00	386	2,38	
		4,00	686	2,97	
		8,00	1449	2,53	
<i>FLC lambda</i>			Date	10/04/06	
Reagent	R.IMG.FRL	Lot	IMFL-A603		
Calibrator	W.UPC.FRL	Lot	WFL-A603		
	Procedure	Standard	HS		
	mg/dl	signal	CV%	signal	CV%
	Saline	0	110	2,02	
		0,125	126	1,11	
		0,250	153	2,19	
		0,50	179	2,79	
		1,00	259	0,77	
		2,00	449	1,06	
		4,00	836	1,28	
		8,00	1525	0,75	