Italiano

Scopo

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English

Indented Use

The Retinol Binding Protein (RBP) "IMG" Kits, when used in conjunction with Beckman Coulter Immage® Analyzer, are intended for the determination of Retinol Binding Protein in serum and urine.

Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with antibodies specific to RBP.

Kit Type		serum/urine	KR	K
<i>Antiserum Reagents</i>				
<input type="checkbox"/> Retinol Binding Protein serum	<i>R.IMG.RBP-S</i>	1 vial – 3,5 ml	serum	serum
<input type="checkbox"/> Retinol Binding Protein urine	<i>R.IMG.RBP-U</i>	1 vial – 4,7 ml	urine	urine
<i>Accessory Reagent "B"</i>				
<input type="checkbox"/> Reagent "B" RBP	<i>R.IMG.RBP.B</i>	1 vial – 1 ml	•	•
<i>Calibrators</i>				
<input type="checkbox"/> Retinol Binding Protein	<i>W.SPC.RBP</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Prediluted Calibrator Set RBP <p>A Set of ready-to-use Prediluted Calibrators is available.</p>	<i>W.SET.SPC.RBP</i>	5 vials – 1 ml	-	-
<i>Controls</i>				
<input type="checkbox"/> Retinol Binding Protein Control Medium Value – Serum	<i>Y.SPR.RBP.MV</i>	1 vial – 1 ml	-	serum
<input type="checkbox"/> Retinol Binding Protein Control Low Value – Serum	<i>Y.SPR.RBP.LV</i>	1 vial – 1 ml	-	serum
<input type="checkbox"/> Retinol Binding Protein Control Medium Value – Urine	<i>Y.UPR.RBP.MV</i>	1 vial – 1 ml	-	urine
<input type="checkbox"/> Retinol Binding Protein Control Low Value - Urine	<i>Y.UPR.RBP.LV</i>	1 vial – 1 ml	-	urine
<ul style="list-style-type: none">Reagents included in the kit. Correlated reagents not included in the kit. They are available separately.				

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

Information re Calibrators / Controls

Due to lack of previsions from the Scientific Committees our "Internal Standardisation" has been applied.

See the relevant leaflet enclosed.

Analysr Set up

Follow the manufacturer's instructions to program one "UDR" assay. The assay parameters can be found on the back.

Getting the cartridges ready

One UDR cartridge must be labelled and filled as above:

- Compartment A: with Antiserum Reagents (*R.IMG.RBP.S/U*).
- Compartment B: with Accessory Reagent "B" (*R.IMG.RBP.B*).

Processing Procedure

Calibration Curve

Prepare various calibrator levels (see table); calibrators dilutions are stable for a few hours.

As an alternative use the Prediluted Calibrator Set – W.SET.SPC.RBP.

Calibration Procedure

The calibration procedure must be run for the chemistry.

The calibration curve must be rerun at every cartridge replacement and whenever a new lot number is used

Results and curves are reported on the back.

Sample

- Serum: Fresh serum, limpид, not hemolyzed, not lipemic.
- Urine: Unconcentrated urine, centrifuged (3.000 rpm x 10'). Use 0,1% Sodium Azide as preservative.

Assay on samples

Samples must be run as provided for by analyser instruction.

Quality Control

	serum	urine
<input type="checkbox"/> Control Medium value	Y.SPR.RBP.MV	Y.UPR.RBP.MV
<input type="checkbox"/> Control Low value	Y.SPR.RBP.LV	Y.UPR.RBP.LV

Results Interpretation

Reference Range

It is advisable that each laboratory settles its own normal reference range.

The literature report

- serum: 3 – 6 mg/dl
- urine: < 0,1 mg/dl

⚠️ Advice on Reagent usage

100 tests can be run on one reagent vial; please use serological pipette to avoid reagent leaks filling cartridges.

Cartridges Compartment dead volume are 250 µl.

⚠️ Precautions and Hazards

Preservative: Sodium Azide.

- Cautlon: Sodium Azide is toxic.
 - All reagents contain Sodium Azide < 0,1% (w/v). **Sodium Azide is not dangerous at this concentration.** Careful handling is recommended to avoid accidental swallowing.
- Waste: Sodium Azide may form explosive compounds in metal drainage.

⚠️ Warnings

- All reagents are intended for *in vitro "research"*.
- Reagents have passed Quality Control checks and are expected to react as described above. Should controls samples yield unexpected results, reagents are not to be considered reliable.
- Do not mix reagents from different lot numbers.
- Calibrators dilutions are stable for a few hours.

Technical Information

Please call Customer Support Service or visit www.newsscientific.com.

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