## Antigen Excess

The proposed application allows an accurate assay of a calibrator concentration four times higher than the highest set point of the calibration curve. However, "Antigen Excess" cannot be ruled out.

## ADVICE ON REAGENT CONSUMPTION

The volume of reagents supplied in the vials allows to run 200 test with the proposed standard parameters; to avoid loss filling the reagent containers it is recommended to use a fine point pipette.

## PRECAUTIONS AND HAZARDS

- **Preservative:** Sodium Azide.
  - **Handling:** Sodium Azide is toxic. All the reagents contain Sodium Azide <0,1% (w/v). Even if Sodium Azide is not dangerous at this concentration, make sure to take the necessary precautions to prevent accidental swallowing.
  - **Waste:** Sodium Azide can react with lead and copper plumbing to form highly explosive metal azides.
- Since absence of infectious agents can not be proven with absolute certainty, all materials derived from human or animal blood should always be handled with due care, following the precautions recommended for biohazardous material.

## WARNINGS

- All reagents are intended only for "In Vitro research".
- Reagents have passed quality control checks and are expected to react as described in this Instructions. If controls do not yield the predicted reaction, all the reagents would not be considered reliable.
- Do not mix reagents from different lot kits.

## TECHNICAL INFORMATION

Contact Customer Support Service or visit [www.newscientific.com](http://www.newscientific.com).

## TECHNICAL PARAMETERS

(rev.: 07.01.2004)

### CALIBRATION CURVE DILUTIONS

LEVEL	DILUTION	FACTOR (% NEAT)	VOL. CALIBRATOR	VOL. SALINE
# 1	1 : 8	12.5 %	100 ul	700 ul
# 2	1 : 4	25 %	100 ul	300 ul
# 3	1 : 2	50 %	100 ul	100 ul
# 4	3 : 4	75 %	300 ul	100 ul
# 5	1 : 1	100 %	200 ul	---

### ASSAY PARAMETERS SET UP - GENERAL INFORMATION

#### Option A: Without Blank Reagent (use only Antiserum Reagent)

**Volumes**  
SAMPLE: 10.0 ul (PREDILUTED 1:20)  
REAGENT: 250 ul (ANTISERUM REAGENT)

#### Measurements

WAVELENGTH: 340 nm  
POINT 1 (BLANK): as soon as possible after mixing sample and reagent (from 0 to (maximum) 10 seconds)  
POINT 2 ("FINAL" POINT): from 8 to 10 minutes, after mixing sample and reagent

#### Calibration Parameters

CALIB. TYPE: Not Lineal (5 points)  
FORMULA: 3rd Order Polinomial, Logit or Polygonal

#### Calculation

Interpolate the absorbance increment (P2-P1) of the samples and controls in the curve obtained with the calibrators

#### Option B: With Blank Reagent (use the two Reagents)

**Volumes**  
SAMPLE: 10.0 ul (PREDILUTED 1:20)  
REAGENT : 250 ul (ANTISERUM OR BLANK REAGENT)

#### Measurements

WAVELENGTH: 340 nm  
POINT 1 (BLANK): after dispensing reagent and just before dispensing sample, or as soon as possible after mixing sample and reagent  
POINT 2 ("FINAL" POINT): from 8 to 10 minutes, after mixing sample and reagent

#### Calibration Parameters

CALIB. TYPE: Not Lineal (5 points)  
FORMULA: 3rd Order Polinomial, Logit or Polygonal

#### Calculation

Interpolate the absorbance increment (P2-P1), obtained with the Main Reagent (Antiserum), minus the absorbance increment (P2-P1), obtained with the Blank Reagent (RWA), of the samples and controls in the curve obtained with the calibrators


### EXAMPLE OF CALIBRATION CURVE - "Cobas Mira" Analyzer

#### Results

REF.: R.ITS.KAP.AS  
LOT: 2156  
CALIBRATOR: REF.: W.SPC.IS1  
LOT: 3371

dil. fact %	Kap-LC mg/dl	Ig-Kap mg/dl	signal mAU
12.5	87	300	0.0271
25	174	600	0.0874
50	348	1200	0.1832
75	522	1800	0.2657
100	696	2400	0.3326

#### Graphic

#### Note - Warning

The results of the calibration and the graphic have value only as an example

BIOMEDICINA - PATOLOGIA CLINICA - INFORMATICA

Immuno Turbidimetric System

**REF** K.ITS.KAP

**Cadenas Ligeras Kappa (Ig-K) - Kit ITS**  
**Kappa Light Chains (Ig-K) - Kit ITS**

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**INTENDED USE**

The *Kappa Light Chains (Ig-K) - Kit ITS* Kit, when used with Clinical Chemistry photometric analyzers, is intended for the determination of **Kappa Light Chains (bound and free)** or **Kappa Immunoglobulins (KAP)** in human serum.

**METHOD**

The method is based on the immunoprecipitation reaction in liquid phase (turbidimetry) with antibodies specific to the Kappa Light Chains of Immunoglobulins.

**CONTENTS - KIT COMPONENTS**

**Prediluted Antiserum Reagent:**

- *Light Chains Kappa (B&F)* (R.ITS.KAP.AS) 2 vials – 25+25 ml

**Reagent Without Antiserum:**

- *Light Chains - RWA* (R.ITS.RWA.LC) 2 vials – 25+25 ml

**Form:** liquid, ready for use.

**Stability:** 12 months at +2°C ... +8°C.

**OTHER RELATED REAGENTS, NOT INCLUDED IN THE KIT**

**Calibrators:**

- *Inmunoglobulins - Calibrator* (W.SPC.IS1) 1 vial – 1 ml

**Controls:**

- *Inmunoglobulins - Medium Value* (W.SPR.IS1.MV) 1 vial – 1 ml
- *Inmunoglobulins - Low Value* (W.SPR.IS1.LV) 1 vial – 1 ml

**Information about Calibrators and Controls:**

See the relevant Instructions.

**ANALYZER SET UP**

Follow the manufacturer of the analyzer instructions to program a monorreactive assay. If the Reagent Without Antiserum is used as "blank reagent" to make a dynamic sample blank, follow the manufacturer of the analyzer instructions to program a monorreactive assay with "blank reagent" (in some analyzers it could be necessary to program a second assay with the same parameters as the main assay). The standard assay parameters can be found on the back.

**OPERATING PROCEDURE**

**Reagent Containers Preparation**

Drain carefully the reagent vials in the proper analyzer reagent containers, using as:

- Main Reagent: Prediluted Antiserum Reagent
- Blank Reagent (optional): Reagent Without Antiserum

**Calibration Curve**

Set up the calibrator (W.SPC.IS1) dilutions as shown in the table on the back.

**Calibration Procedure**

Run the calibration procedure preveded for the analyzer. For the calibration curve construction use a "Non Lineal" model (for example 3rd Order Polynomial, Logit or Poligonal type).

**The calibration must be rerun, at least, with each new lot.**

The results and a typical curve, obtained with a "Cobas Mira" analyzer, are reported as an example on the back.

**Samples**

Fresh Serum, limpid, not hemolyzed, not lipemic.

**Predilute the serum 1:20 for its analysis.**

**Assay on Samples**

Samples must be run as preveded for the analyzer.

**Quality Control**

It is advisable to insert in each analytical serie:

- Medium Value Control: W.SPR.IS1.MV
- Low Value Control: W.SPR.IS1.LV

**RESULTS AND INTERPRETATION**

**Reference Values**

It is advisable that each laboratory settles its own normal reference range. KAP reference values specified in literature are between 200-440 mg/dl (expressed as Kappa Light Chain concentration) or between 650-1500 mg/dl (expressed as Kappa Immunoglobulin concentration).

Kappa/Lambda ratio reference values specified in literature are between 1,35-2,65.

**Limitations – Lack of Sample/Calibrator Parallelism**

The Light Chains in the sample could be "different" from those in calibrators and therefore they could react unproportionately; *lack of sample/calibrator parallelism*. The same occurs with the quantification of serum monoclonal immunoglobulins (serum Monoclonal Component – MC) by liquid phase Immunoprecipitation (Turbidimetry, Nephelometry).