Monoclonality

This test does not allow highlighting if the Light Chains in the sample are monoclonal or polyclonal.

ADVICE ON REAGENT CONSUMPTION

The volume of reagents supplied in the vials allows to run 100 test with the proposed standard parameters; to avoid loss filling the reagent containers it is recommended to use a fine point pipette.

PRECAUTIONS AND HAZARDS

- **Preservative:** Sodium Azide.
 - **Handling:** Sodium Azide is toxic.
All the reagents contain Sodium Azide <0,1% (w/v). Even if Sodium Azide is not dangerous at this concentration, make sure to take the necessary precautions to prevent accidental swallowing.
 - **Waste:** Sodium Azide can react with lead and copper plumbing to form highly explosive metal azides.
- Since absence of infectious agents can not be proven with absolute certainty, all materials derived from human or animal blood should always be handled with due care, following the precautions recommended for biohazardous material.

WARNINGS

- All reagents are intended only for *In Vitro* tests.
- Reagents have passed quality control checks and are expected to react as described in this Instructions. If controls do not yield the expected reaction, all the reagents should not be considered reliable.
- Do not mix reagents from different lot kits.
- Calibrator dilutions are stable only for their immediate use.
- With the pass of time or unsuitable conservation conditions, for example with sudden or frequent temperature changes, the reagents can form a slight inespecific precipitate. Filtering with a 0.2-0.45µm monouse syringe filter, the quality and performances of the reagents remain unchanged.

TECHNICAL INFORMATION

Contact Customer Support Service or visit www.newscientific.com.

TECHNICAL PARAMETERS

(rev.: 03.08.2010)

CALIBRATION CURVE DILUTIONS

LEVEL	DILUTION	FACTOR (% NEAT)	VOL. CALIBRATOR	VOL. SALINE
# 1	1 : 16	6.25 %	50 ul	750 ul
# 2	1 : 8	12.5 %	50 ul	350 ul
# 3	1 : 4	25 %	100 ul	300 ul
# 4	1 : 2	50 %	200 ul	200 ul
# 5	1 : 1	100 %	400 ul	---

ASSAY PARAMETERS SET UP - GENERAL INFORMATION

Volumes

SAMPLE/CALS./CTS.: 10.0 ul (PRE-DILUTED 1:15)
REAGENT-1: 250 ul (BUFFER - R.ITS.TLC.BF)
REAGENT-2: 50 ul (ANTISERUM REAGENT - STARTER - R.ITS.LAM.AS)

Measurements

PRIMARY WAVELENGTH: 340 nm
SECONDARY WAVELENGTH: NONE OR >600 nm
POINT 1 (BLANK): a few seconds before dispensing R2 (as starter reagent) and at least 2-3 minutes after dispensing and mixing sample and R1 (buffer).
POINT 2 ("FINAL" POINT): at a fixed time (between 3 to 10 minutes), after dispensing and mixing R2

Calibration Parameters

CALIB. TYPE: Not Lineal (5 points)
FORMULA: 3rd Order Polynomial, Logit or Polygonal (select the best adjustment)

Calculation

Interpolate the absorbance increment (P2-P1) of the samples and controls in the curve obtained with the calibrators


EXAMPLE OF CALIBRATION CURVE - "Cobas Mira" Analyzer

Results

REAGENTS: REF.: KR.ITS.LAM
LOT: J081
CALIBRATOR: REF.: W.SPC.IS1
LOT: J081

dil. fact. %	Lam-LC mg/dl	Ig-Lam mg/dl	signal mAU
6.25	40.8	140	0,04620
12.5	81.6	280	0,05935
25	163	559	0,10665
50	326	1118	0,22245
100	653	2237	0,43720

Graphic

Note - Warning

The results of the calibration and the graphic have value only as an example

BIOMEDICINA - PATOLOGIA CLINICA - INFORMATICA

Immuno Turbidimetric System

KR.ITS.LAM

Cadenas Ligeras Lambda (Ig-L) - Kit Reactivos ITS
Lambda Light Chains (Ig-L) - ITS Reagent Kit

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INTENDED USE

The "*Lambda Light Chains (Ig-L) - ITS Reagent Kit*" Kit, when used with Clinical Chemistry photometric analyzers, is intended for the determination of **Lambda Light Chains (bound and free) or Lambda Immunoglobulins (LAM)** in human serum.

METHOD

The method is based on the immunoprecipitation reaction in liquid phase (turbidimetry/nephelometry) with antibodies specific to the Lambda (B&F) Light Chains of Immunoglobulins.

CONTENTS - KIT COMPONENTS

Antiserum Reagent:

- Light Chains Lambda (B&F) - Antiserum Reagent* (R.ITS.LAM.AS) 1 vial – 5 ml

Reaction Buffer:

- Light Chains (B&F) - Buffer* (R.ITS.TLC.BF) 1 vial – 25 ml

Form: liquid, ready for use.

Stability: 12 months at +2°C ... +8°C.

OTHER RELATED REAGENTS, NOT INCLUDED IN THE KIT

Calibrators:

- Inmunoglobulins - Calibrator* (W.SPC.IS1) 1 vial – 1 ml

Controls:

- Inmunoglobulins - Medium Value* (Y.SPR.IS1.MV) 1 vial – 1 ml
- Inmunoglobulins - Low Value* (Y.SPR.IS1.LV) 1 vial – 1 ml

Information about Calibrators and Controls:

See the relevant Instructions.

ANALYZER SET UP

Follow the manufacturer of the analyzer instructions to program a birreactive assay.

The standard assay parameters can be found on the back.

OPERATING PROCEDURE

Reagent Containers Preparation

Drain carefully the reagent vials in the proper analyzer reagent containers, using as:

- Reagent 1 (buffer): Reaction Buffer (R.ITS.TLC.BF)
- Reagent 2 (starter): Antiserum Reagent (R.ITS.LAM.AS)

Calibration Curve

Set up the calibrator (W.SPC.IS1) dilutions as shown in the table on the back.

Calibration Procedure

Run the calibration procedure foreseen for the analyzer.

For the calibration curve construction use a "Non Linear" model (for example 3rd Order Polynomial, Logit or Polygonal type).

The calibration must be rerun, at least, with each new lot.

The results and a typical curve, obtained with a "Cobas Mira" analyzer, are reported as an example on the back.

Samples

Fresh Serum, límpid, not hemolyzed, not lipemic.

Predilute the serum 1:15 for its analysis.

Assay on Samples

Samples must be run as foreseen for the analyzer.

Quality Control

It is advisable to insert in each analytical serie:

- Medium Value Control: Y.SPR.IS1.MV
- Low Value Control: Y.SPR.IS1.LV

RESULTS AND INTERPRETATION

Reference Values

It is advisable that each laboratory settles its own normal reference range.

LAM reference values specified in literature are between 110-240 mg/dl (expressed as Lambda Light Chain concentration) or between 350-800 mg/dl (expressed as Lambda Immunoglobulin concentration).

Kappa/Lambda ratio reference values specified in literature are between 1,35-2,65.

Limitations – Lack of Sample/Calibrator Parallelism

The Light Chains in the sample could be "different" from those in calibrators and therefore they could react unproportionately; *lack of sample/calibrator parallelism*. The same occurs with the quantification of serum monoclonal immunoglobulins (serum Monoclonal Component – MC) by liquid phase Immunoprecipitation (turbidimetry/nephelometry).

Antigen Excess

The proposed application allows an accurate assay of a calibrator concentration four times higher than the highest set point of the calibration curve. However, "Antigen Excess" cannot be ruled out.