PRECAUTIONS AND HAZARDS

- **Preservative:** Sodium Azide.
- **Handling:** Sodium Azide is toxic. All the reagents contain Sodium Azide <0,1% (w/v). Even if Sodium Azide is not dangerous at this concentration, make sure to take the necessary precautions to prevent accidental swallowing.
- **Waste:** Sodium Azide can react with lead and copper plumbing to form highly explosive metal azides.
- Since absence of infectious agents can not be proven with absolute certainty, all materials derived from human or animal blood should always be handled with due care, following the precautions recommended for biohazardous material.

WARNINGS

- All reagents are intended only for *In Vitro* assays.
- Reagents have passed quality control checks and are expected to react as described in this instructions. If controls do not yield the preveved reaction, all the reagents would not be considered reliable.
- Do not mix reagents from different lot kits.

TECHNICAL INFORMATION

Contact Customer Support Service or visit www.newscientific.com.

TECHNICAL PARAMETERS

(rev.: 15.03.2005)

CALIBRATION CURVE DILUTIONS

LEVEL	DILUTION	FACTOR (% NEAT)	VOL. CALIBRATOR	VOL. SALINE
# 1	- - -	0 %	- - -	300 ul
# 2	1 : 12	8.33 %	100 ul	1100 ul
# 3	1 : 6	16.66 %	100 ul	500 ul
# 4	1 : 3	33.33 %	100 ul	200 ul
# 5	2 : 3	66.66 %	200 ul	100 ul
# 6	1 : 1	100 %	300 ul	- - -

ASSAY PARAMETERS SET UP - GENERAL INFORMATION

Volumes

SAMPLE (NOT PREDILUTED): 2.0 ul
R1 (REACTION BUFFER): 250 ul
R2 (LATEX AS. REAGENT): 60 ul

Measurements

WAVELENGTH: 500 nm
POINT 1 (SAMPLE BLANK): after mixing Sample and Reagent 1 (reaction buffer)
POINT 2 ("FINAL" POINT): from 5 to 10 minutes, after dispensing Reagent 2 (latex as. reagent (start reagent))

Calibration Parameters

CALIB. TYPE: Not Lineal (6 points)
FORMULA: 3rd Order Polynomial, Logit or Polygonal

Procedure - Calculation

- Mix the Sample and the Reagent 1
- Incubate a fixed time (at least 1 minute)
- Read the absorbance (P1)
- Dispense the Reagent 2
- Incubate a fixed time (between 5 and 10 minutes)
- Read the absorbance (P2)
- Interpolate the absorbance increment (P2-P1) of the samples and controls in the curve obtained with the calibrators

BIOMEDICINA - PATOLOGIA CLINICA - INFORMATICA

Immuno Turbidimetric System (Latex)

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INTENDED USE

The *Alpha-1 Microglobulin - urine - Kit ITSL* Kit, when used with Clinical Chemistry photometric analyzers, is intended for the determiniarion of **Alpha-1 Microglobulin (A1M)** in human urine.

METHOD

The method is based on the immunoprecipitation reaction in liquid phase (turbidimetry) with antibodies specific to A1M, bound with latex particles.

CONTENTS - KIT COMPONENTS

Latex Antiserum Reagent:

- Alpha-1 Microglobulin - Latex* (R.ITSLA1M.AS) 1 vial – 6 ml

Reaction Buffer:

- Alpha-1 Microglobulin - Buffer* (R.ITSLA1M.BF) 1 vial – 25 ml

Form: liquid, ready for use.

Stability: 12 months at +2°C ... +8°C.

OTHER RELATED REAGENTS, NOT INCLUDED IN THE KIT

Calibrator:

- Alpha-1 Microglobulin - Calibrator* (W.A1M.STD) 1 vial – 1.5 ml

Control:

- Alpha-1 Microglobulin - Control* (W.A1M.CT) 3 vials – 1 ml

Information about Calibrators and Controls:

See the relevant Instructions.

ANALYZER SET UP

Follow the manufacturer of the analyzer instructions to program a birreactive assay.

The standard assay parameters can be found on the back.

OPERATING PROCEDURE

Reagent Containers Preparation

Drain carefully the reagent vials in the proper analyzer reagent containers, using as:

- Reagent 1: Buffer Reaction
- Reagent 2: Latex Antiserum Reagent

Calibration Curve

Set up the calibrator (W.A1M.STD) dilutions as shown in the table on the back.

Calibration Procedure

Run the calibration procedure preveded for the analyzer.

For the calibration curve construction use a "Non Lineal" model (for example 3rd Order Polynomial, Logit or Poligonal type).

The calibration must be rerun, at least, with each new lot.

Samples

Fresh Urine, limpid, centrifuged (3000 rpm - 10')..

Assay on Samples

Samples must be run as preveded for the analyzer.

Quality Control

It is advisable to insert in each analittical serie the control - *W.A1M.CT*.

RESULTS AND INTERPRETATION

Reference Values

It is advisable that each laboratory settles its own normal reference range.

A1M reference values specified in literature usually are expressed in relation to the diuresis (mg/24h, µg/min or refered to the Creatinine excretion). An excretion higher than 14 mg/24h is usually considered pathologic.

ADVICE ON REAGENT CONSUMPTION

The volume of reagents supplied in the vials allows to run 100 test with the proposed standard parameters; to avoid loss filling the reagent containers it is recommended to use a fine point pipette.