ADVICE ON REAGENT CONSUMPTION

The volume of reagents supplied in the vials allows to run 100 test with the proposed standard parameters; to avoid loss filling the reagent containers it is recommended to use a fine point pipette.

PRECAUTIONS AND HAZARDS

- **Preservative:** Sodium Azide.
 - **Handling:** Sodium Azide is toxic. All the reagents contain Sodium Azide <0,1% (w/v). Even if Sodium Azide is not dangerous at this concentration, make sure to take the necessary precautions to prevent accidental swallowing.
 - **Waste:** Sodium Azide can react with lead and copper plumbing to form highly explosive metal azides.
- Since absence of infectious agents can not be proven with absolute certainty, all materials derived from human or animal blood should always be handled with due care, following the precautions recommended for biohazardous material.

WARNINGS

- All reagents are intended only for *In Vitro* tests.
- Reagents have passed quality control checks and are expected to react as described in this Instructions. If controls do not yield the preveved reaction, all the reagents would not be considered reliable.
- Do not mix reagents from different lot kits.

TECHNICAL INFORMATION

Contact Customer Support Service or visit www.newscientific.com.

ASSAY PARAMETERS SET UP - GENERAL INFORMATION

Volumes

SAMPLE (NOT PREDILUTED): 7.0 ul
 R1 (REACTION BUFFER): 200 ul
 R2 (LATEX A. REAGENT): 75 ul

Measurements

WAVELENGTH: 600 nm
 POINT 1 (BLANK): around 30 seconds after mixing sample and reagents (usually, a measurement between 20 and 40 seconds after mixing is better than an immediate measurement because the mixture is more stabilized). As an alternative, the 1st measurement (blank) can be done few seconds before dispensing R2 (as the starter reagent), 5 minutes, after mixing sample and reagents

POINT 2 ("FINAL" POINT):

Calibration Parameters


CALIB. TYPE: Not Lineal (6 points)
 FORMULA: 3rd Order Polinomial, Logit or Polygonal

Calculation

Interpolate the absorbance increment (P2-P1) of the samples and controls in the curve obtained with the calibrators

EXAMPLE OF CALIBRATION CURVE - "Cobas Mira" Analyzer

REAGENTS: REF.: K.ITSL.IGE
 LOT: I021
 CALIBRATOR: REF.: W.IGE.STD-K
 LOT: H051

Note - Warning

The results of the calibration and the graphic have value only as an example

BIOMEDICINA - PATOLOGIA CLINICA - INFORMATICA

Immuno Turbidimetric System (Latex)

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METHOD

The method is based on the immunoprecipitation reaction in liquid phase (turbidimetry) with antibodies specific to IGE, bound with latex particles.

CONTENTS - KIT COMPONENTS

Latex Antiserum Reagent:

- IgE Inmunoglobulins - Latex* (R.ITSL.IGE.AS) 1 vial – 7.5 ml

Reaction Buffer:

- IgE Inmunoglobulins - Buffer* (R.ITSL.IGE.BF) 1 vial – 20 ml

Form: liquid, ready for use.

Stability: more than 12 months at +2°C ... +8°C.

OTHER RELATED REAGENTS, NOT INCLUDED IN THE KIT

Calibrators:

- IgE Inmunoglobulins - Calibrator Set* (W.IGE.STD-K) 6 vials (6 levels) – 1 ml

Controls:

- IgE Inmunoglobulins - Control Set* (W.IGE.CT-K) 6 vials (2 levels) – 1 ml

Information about Calibrators and Controls:

See the relevant Instructions.

ANALYZER SET UP

Follow the manufacturer of the analyzer instructions to program a birreactive assay. The standard assay parameters can be found on the back.

OPERATING PROCEDURE

Reagent Containers Preparation

Drain carefully the reagent vials in the proper analyzer reagent containers, using as:

- Reagent 1: Buffer Reaction
- Reagent 2: Latex Antiserum Reagent

Calibration Curve

The use of the six levels prediluted calibrators set – *W.IGE.STD-K*, is recommended.

Calibration Procedure

Run the calibration procedure preveded for the analyzer.

For the calibration curve construction use a "Non Lineal" model (for example 3rd Order Polinomial, Logit or Poligonal type).

The calibration must be rerun, at least, with each new lot.

The results and a tipycal curve, obtained with a "Cobas Mira" analyzer, are reported as an example on the back.

Samples

Fresh Serum, limpid, not hemolyzed, not lipemic.

Assay on Samples

Samples must be run as preveded for the analyzer.

Quality Control

It is advisable to insert in each analytical serie the controls of the two level control set - *W.IGE.CT-K*.

RESULTS AND INTERPRETATION

Reference Values

It is advisable that each laboratory settles its own normal reference range.

The literature reports the following values:

- newborn: up to 1.5 IU/ml
- children: up to 1 year: up to 15 IU/ml
 - 1 - 5 years: up to 60 IU/ml
 - 6 - 9 years: up to 90 IU/ml
 - 10 - 15 years: up to 200 IU/ml
- adults: up to 100 IU/ml