Italiano

Scopo

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English

Indented Use

The C1q (C1Q) "APS" Kits, when used in conjunction with Beckman Coulter "APS-Array [®]" Analyzer, are intended for the determination of C1q in serum

Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with antibodies specific to C1q.

Reagents

Kit Type			KR	K
<i>Antiserum Reagent</i>				
<input type="checkbox"/> C1 (esterase) inhibitor	<i>R.APS.C1Q</i>	1 vial – 4,5 ml	•	•
<i>Calibrator</i>				
<input type="checkbox"/> Complement system protein	<i>W.SPC.CS1</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Prediluted Calibrator Set Complement system proteins A Set of ready-to-use Prediluted Calibrators is available.	<i>W.SET.SPC.CS1</i>	5 vials – 1 ml	-	-
<i>Controls</i>				
<input type="checkbox"/> Complement system proteins Medium Value	<i>Y.SPR.CS1.MV</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Complement system proteins Low Value	<i>Y.SPR.CS1.LV</i>	1 vial – 1 ml	-	•
• Reagents included in the kit. - Correlated reagents not included in the kit. They are available separately.				

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

Information re Calibrators / Controls

Due to lack of previsions from the Scientific Committees our "Internal Standardisation" has been applied.

See the relevant leaflet enclosed.

Analysér Set up

Follow the manufacturer's instructions to program a "UDR" assay. The assay parameters can be found on the back.

Processing Procedure

Calibration Curve

Prepare various calibrator levels (see table): calibrators dilutions are stable for a few hours.

As an alternative use the Prediluted Calibrator Set – W.SET.SPC.CS1.

Calibration Procedure

The calibration procedure must be run for the chemistry.

The calibration curve must be rerun whenever a new lot number is used.

Put the concentrations of the calibration points on axis-x and the corresponding "Results" (Reagent Blank net) on axis-y. Results and curves are reported on the back.

Sample

- Serum: Fresh serum, limpid, not hemolyzed, not lipemic.

Assay on samples

Samples must be run as provided for by analyser instruction.

To calculate sample concentration interpolate sample Rate result (Reagent Blank net) on the Calibration curve.

Quality Control

For each set of determinations it is recommended to include:

<input type="checkbox"/> Control Medium value	Y.SPR.CS1.MV
<input type="checkbox"/> Control Low value	Y.SPR.CS1.LV

Results Interpretation

Reference Range

It is advisable that each laboratory settles its own normal reference range.

The literature report: 10 – 25 mg/dl (Putnam F.W. – The Plasma Protein – II Edition).

⚠️ Advice on Reagent usage

100 tests can be run on one reagent vial. Reagent vial dead volume is 400 µl.

⚠️ Precautions and Hazards

Preservative: SodioAzide.

- Caution: Sodium Azide is toxic. All reagents contain Sodium Azide < 0,1% (w/v). Sodium Azide is not dangerous at this concentration. Careful handling is recommended to avoid accidental swallowing.
- Waste: Sodium Azide may form explosive compounds in metal drainage.

⚠️ Warnings

- All reagents are intended for *in vitro "research"*.
- Reagents have passed Quality Control checks and are expected to react as described above. Shold controls samples yield unexpected results, reagents are not to be considered reliable.
- Do not mix reagents from different lot numbers.
- Calibrators dilutions are stable for a few hours.

Technical Information

Please call Customer Support Service or visit www.newscientific.com.

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Technical Parameters

C1Q - Calibration Curve Dilution


REF	W.SPC.C1Q	LOT	A311				
Prepare various calibrator levels doubling or as shown in the table below:							
Calibrator Dilutions			1 : 16	1 : 8	1 : 4	1 : 2	1 : 1
<input type="checkbox"/>	Calibrator	µl	25	50	50	100	200
<input type="checkbox"/>	Buffered Saline (PBS)	µl	375	350	150	100	0
Equivalent Calibrator Concentration - Neat mg/dl (*)		30	1,88	3,75	7,50	15,00	30
* The calibrator values are lot dependent.							
As an alternative use Prediluted Calibrator Set: Code W.SET.SPC.CS1							

Assay Parameters set up on APS

Program the assay parameters according to the following			
UPC ID	C1Q		
Reagent position	Select free position		
Gain	44		
Optics	L (left) or R (right)		
Dilution	Calibration curve	1	: 6
	Samples	1	: 6

C1Q - Calibration Curve Example - APS

<i>C1Q</i>		Date	03/11/2003	
Reagent	R.APS.C1Q	Lot	A310	
Calibrator	W.SPC.SC2	Lot	N14801-00	
Net	mg/dl	30		
pre Dilution	mg/dl	signal	CV%	Sample mg/dl
1 : 16	1,88	70		1,88
1 : 8	3,75	177		3,75
1 : 4	7,50	429		7,50
1 : 2	15,00	1220		15,00
1 : 1	30,00	2360		30,00

Please Note: the value of the curves shown are only indicative

BIOMEDICINA – PATOLOGIA CLINICA - INFORMATICA

ImmunoNefelometria su APS-Array®
ImmunoNephelometry on APS-Array®
ImmunNephelometrie auf APS-Array®

K.APS.C1Q
KR.APS.C1Q


Reagenti per uso professionale
Reagents for professional use
Reagenzien für Berufsgebrauch

Complemento C1q
Complement C1q
Complement C1q

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New Scientific Company S.r.l. – Via Dante Alighieri, 35 – I 20032 Cormano (MI) – Italy