# Italiano

## Scopo

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# English

## Indented Use

The C1 (esterase) inhibitor (CEI) "APS" Kits, when used in conjunction with Beckman Coulter "APS-Array <sup>®</sup>" Analyzer, are intended for the determination of C1 (esterase) inhibitor in serum.

## Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with antibodies specific to CEI.

## Reagents

Kit Type			KR	K
<b>Antisera Reagent</b>				
<input type="checkbox"/> C1 (esterase) inhibitor	<i>R.APS.CEI</i>	1 vial – 4,5 ml	•	•
<b>Calibrator</b>				
<input type="checkbox"/> Complement system protein	<i>W.SPC.CS1</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Prediluted Calibrator Set Complement system proteins A Set of ready-to-use Prediluted Calibrators is available.	<i>W.SET.SPC.CS1</i>	5 vials – 1 ml	-	-
<b>Controls</b>				
<input type="checkbox"/> Complement system proteins Medium Value	<i>Y.SPR.CS1.MV</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Complement system proteins Low Value	<i>Y.SPR.CS1.LV</i>	1 vial – 1 ml	-	•
• Reagents included in the kit. - Correlated reagents not included in the kit. They are available separately.				

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

## Information re Calibrators / Controls

Due to lack of provisions from the Scientific Committees our "Internal Standardisation" has been applied.

See the relevant leaflet enclosed.

## Analysers Set up

Follow the manufacturer's instructions to program a "UDR" assay. The assay parameters can be found on the back.

## Processing Procedure

### Calibration Curve

Prepare various calibrator levels (see table): calibrators dilutions are stable for a few hours.

As an alternative use the Prediluted Calibrator Set – W.SET.SPC.CS1.

### Calibration Procedure

The calibration procedure must be run for the chemistry.

The calibration curve must be rerun whenever a new lot number is used.

Put the concentrations of the calibration points on axis-x and the corresponding "Results" (Reagent Blank net) on axis-y. Results and curves are reported on the back.

### Sample

Serum: Fresh serum, limpid, not hemolyzed, not lipemic.

### Assay on samples

Samples must be run as provided for by analyser instruction.

To calculate sample concentration interpolate sample rate result (Reagent Blank net) on the Calibration curve.

### Quality Control

For each set of determinations it is recommended to include:

<input type="checkbox"/> Control Medium value	Y.SPR.CS1.MV
<input type="checkbox"/> Control Low value	Y.SPR.CS1.LV

## Results Interpretation

### Reference Range

It is advisable that each laboratory settles its own normal reference range.

The literature report: 15 – 35 mg/dl (Putnam F.W. – The Plasma Protein – II Edition).

## ⚠ Advice on Reagent usage

100 tests can be run on one reagent vial. Reagent vial dead volume is 400 µl.

## ⚠ Precautions and Hazards

### Preservative: Sodium Azide.

• Caution: Sodium Azide is toxic.

All reagents contain Sodium Azide < 0,1% (w/v). Sodium Azide is not dangerous at this concentration. Careful handling is recommended to avoid accidental swallowing.

• Waste: Sodium Azide may form explosive compounds in metal drainage.

## ⚠ Warnings

• All reagents are intended for *in vitro* "research".

• Reagents have passed Quality Control checks and are expected to react as described above.

Should controls samples yield unexpected results, reagents are not to be considered reliable.

• Do not mix reagents from different lot numbers.

• Calibrators dilutions are stable for a few hours.

## Technical Information

Please call Customer Support Service or visit [www.newscientific.com](http://www.newscientific.com).

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
## Technical Parameters - APS

CEI - Calibration Curve Dilution - APS					
<b>REF</b>	W.SPC.CS1	<b>LOT</b>	A311		
Prepare various calibrator levels doubling or as shown in the table below:					
Calibrator Dilutions	1 : 16	1 : 8	1 : 4	1 : 2	1 : 1
<input type="checkbox"/> Calibrator	µl	25	50	100	200
<input type="checkbox"/> Buffered Saline (PBS)	µl	375	350	150	0
Equivalent Calibrator					
Concentration - Neat mg/dl (*)	66,6	4,16	8,33	16,65	33,30
* The calibrator values are lot dependent. As an alternative use Prediluted Calibrator Set: Code W.SET.SPC.CS1					


Assay Parameters set up on APS			
Program the assay parameters according to the following			
UPC ID	CEI		
Reagent position	Select free position		
Gain	22		
Optics	L (left) or R (right)		
Dilution	Calibration curve	1	6
	Samples	1	6

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New Scientific Company S.r.l. – Via Dante Alighieri, 35 – I 20032 Cormano (MI) – Italy