# Italiano

## Scopo

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# English

## Indented Use

The Retinol Binding Protein (RBP) “APS” Kits, when used in conjunction with Beckman Coulter “APS-Array ®” Analyzer, are intended for the determination of Retinol Binding Protein in serum and urine.

## Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with antibodies specific to RBP.

## Reagents

Kit Type		serum/urine	KR	K
<b>Antiserum Reagent</b>				
<input type="checkbox"/> Retinol Binding Protein serum / urine	<i>R.APS.RBP</i>	1 vial – 4,5 ml	• serum • urine	• serum • urine
<b>Calibrator</b>				
<input type="checkbox"/> Retinol Binding Protein	<i>W.SPC.RBP</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Prediluted Calibrator Set RBP	<i>W.SET.SPC.RBP</i>	5 vials – 1 ml	A	-
A Set of ready-to-use Prediluted Calibrators is available.				
<b>Controls</b>				
<input type="checkbox"/> Retinol Binding Protein Control Medium Value – Serum	<i>Y.SPR.RBP.MV</i>	1 vial – 1 ml	-	• serum
<input type="checkbox"/> Retinol Binding Protein Control Low Value – Serum	<i>Y.SPR.RBP.LV</i>	1 vial – 1 ml	-	• serum
<input type="checkbox"/> Retinol Binding Protein Control Medium Value – Urine	<i>Y.UPR.RBP.MV</i>	1 vial – 1 ml	-	• urine
<input type="checkbox"/> Retinol Binding Protein Control Low Value - Urine	<i>Y.UPR.RBP.LV</i>	1 vial – 1 ml	-	• urine
• Reagents included in the kit.				
- Correlated reagents not included in the kit. They are available separately.				

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

## Information re Calibrators / Controls

Due to lack of previsions from the Scientific Committees our “Internal Standardisation” has been applied.

See the relevant leaflet enclosed.

## Analysers Set up

Follow the manufacturer's instructions to program a “UDR” assay. The assay parameters can be found on the back.

## Processing Procedure

### Calibration Curve

Prepare various calibrator levels (see table): calibrators dilutions are stable for a few hours.

As an alternative use the Prediluted Calibrator Set – W.SET.SPC.RBP.

### Calibration Procedure

The calibration procedure must be run for the chemistry.

The calibration curve must be rerun whenever a new lot number is used.

Put the concentrations of the calibration points on axis-x and the corresponding “Results” (Reagent Blank net) on axis-y. Results and curves are reported on the back.

### Sample

Serum: Fresh serum, limpid, not hemolyzed, not lipemic.

Urine: Unconcentrated urine, centrifuged (3,000 rpm x 10'). Use 0,1% Sodium Azide as preservative.

### Assay on samples

Samples must be run as provided for by analyser instruction.

To calculate sample concentration interpolate sample Rate result (Reagent Blank net) on the Calibration curve.

### Quality Control

For each set of determinations it is recommended to include:

	serum	urine
<input type="checkbox"/> Control Medium value	Y.SPR.RBP.MV	Y.UPR.RBP.MV
<input type="checkbox"/> Control Low value	Y.SPR.RBP.LV	Y.UPR.RBP.LV

## Results Interpretation

### Reference Range

It is advisable that each laboratory settles its own normal reference range.

The literature report

serum: 3 – 6 mg/dl

urine: < 0,1 mg/dl

## ⚠️ Advice on Reagent usage

100 tests can be run on one reagent vial. Reagent vial dead volume is 400 µl.

## ⚠️ Precautions and Hazards

Preservative: Sodium Azide.

• Caution: Sodium Azide is toxic.

All reagents contain Sodium Azide < 0,1% (w/v). Sodium Azide is not dangerous at this concentration.

Careful handling is recommended to avoid accidental swallowing.

• Waste: Sodium Azide may form explosive compounds in metal drainage.

## ⚠️ Warnings

• All reagents are intended for *in vitro* “research”.

• Reagents have passed Quality Control checks and are expected to react as described above.

Should controls samples yield unexpected results, reagents are not to be considered reliable.

• Do not mix reagents from different lot numbers.

• Calibrators dilutions are stable for a few hours.

## Technical Information

Please call Customer Support Service or visit [www.newscientific.com](http://www.newscientific.com).

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Technical Parameters - Serum and Urine - APS						
RBP - Calibration Curve Dilution - APS						
REF	W.S.P.C.RBP	LOT	A311			
Prepare various calibrator levels doubling or as shown in the table below :						
Calibrator Dilutions		1 : 16	1 : 8	1 : 4	1 : 2	1 : 1
☐ Calibrator	µl	25	50	50	100	200
☐ Buffered Saline (PBS)	µl	375	350	150	100	0
Equivalent Calibrator Concentration - Neat mg/dl (*)		8,02	0,50	1,00	2,01	4,01
* The calibrator values are lot dependent.						
As an alternative use Prediluted Calibrator Set: Code W.SET.SPC.RBP						

Assay Parameters set up on APS			
Program the assay parameters according to the following			
UPC ID	RBS (serum) - RBU (urine)		
Reagent position	Select free position		
Gain	33		
Optics	L (left) or R (right)		
Serum			
Dilution	Calibration curve	1	: 6
	Samples	1	: 6
Urine			
Dilution	Calibration curve	1	: 6
	Samples	1	: 1

RBP - Calibration Curve Example - serum and urine - APS														
RBP - Serum					RBP - Urine									
Reagent	R.APS.RBP	Lot	A311		Date	03/11/2003								
Calibrator	W.S.P.C.RBP	Lot	A311		Reagent	R.APS.RBP	Lot	A311		Date	03/11/2003			
Net	mg/dl	8,02			Net	mg/dl	8,02			Net	mg/dl	8,02		
pre Dilution	mg/dl	signal	CV%	Sample mg/dl	pre Dilution	mg/dl	signal	CV%	Sample mg/dl	pre Dilution	mg/dl	signal	CV%	Sample mg/dl
1 : 16	0,50	62		0,50	1 : 16	0,50	65		0,08	1 : 16	0,50	109		0,17
1 : 8	1,00	108		1,00	1 : 8	1,00	109		0,17	1 : 8	1,00	416		0,33
1 : 4	2,01	407		2,01	1 : 4	2,01	416		0,33	1 : 4	2,01	1440		0,67
1 : 2	4,01	1360		4,01	1 : 2	4,01	1440		0,67	1 : 2	4,01	3620		1,34
1 : 1	8,02	3590		8,02	1 : 1	8,02	3620		1,34	1 : 1	8,02			

Please Note: the value of the curves shown are only indicative

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New Scientific Company S.r.l. - Via Dante Alighieri, 35 - I 20032 Corno (MI) - Italy