Italiano

Scopo

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English

Indented Use

The C1q (C1Q) "BNA" kit , when used in conjunction with Dade Behring "BN™ series" Analyzers (BNA, BN100, BNII, ProSpec), is intended for the determination of C1q in serum.

Method

The method applied is based on ImmunoPrecipitation reaction in liquid phase with antibodies specific to C1q.

Reagents

Kit Type			KR	K
Antiserum Reagent				
<input type="checkbox"/> Complement C1q	<i>R.BNA.C1Q</i>	1 vial – 2 ml	•	•
Calibrator				
<input type="checkbox"/> Complement system protein	<i>W.SPC.CS1</i>	1 vial – 1 ml	-	•
Controls				
<input type="checkbox"/> Complement system proteins Medium Value	<i>Y.SPR.CS1.MV</i>	1 vial – 1 ml	-	•
<input type="checkbox"/> Complement system proteins Low Value	<i>Y.SPR.CS1.LV</i>	1 vial – 1 ml	-	•
• Reagents included in the kit.				
- Correlated reagents not included in the kit. They are available separately.				

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

Information re Calibrators / Controls

Due to lack of provisions from the Scientific Committees our "Internal Standardisation" has been applied.

See the relevant leaflet enclosed.

Analysers Set up

- BNA and BN100** - Follow the manufacturer's instructions to program one new assay. The parameters are on the back.
- BNII and ProSpec** - The protocol is analogous to that for "BNA". For the parameters, contact Customer Services at NSC. For implementation, contact Dade Behring.

Processing Procedure

Calibration Curve

The calibration procedure must be run for each chemistry.

The calibration curve must be rerun whenever a new lot number is used.

Results and curves are reported on the back.

Sample

Serum: Fresh serum, limpid, not hemolyzed, not lipemic.

Assay on samples

Samples must be run as provided for by analyser instruction.

Quality Control

For each set of determinations it is recommended to include:

<input type="checkbox"/> Control Medium value	Y.SPR.CS1.MV
<input type="checkbox"/> Control Low value	Y.SPR.CS1.LV

Results Interpretation

Reference Range

It is advisable that each laboratory settles its own normal reference range.

The literature report: 10 – 25 mg/dl (Putnam F.W. – The Plasma Protein – II Edition).

⚠ Advice on Reagent usage

50 tests can be run on one reagent vial. Reagent vial dead volume is 100 µl.

⚠ Precautions and Hazards

Preservative: Sodium Azide.

- Cautio: Sodium Azide is toxic.
 - All reagents contain Sodium Azide < 0,1% (w/v). **Sodium Azide is not dangerous at this concentration.** Careful handling is recommended to avoid accidental swallowing.
- Waste: Sodium Azide may form explosive compounds in metal drainage.

⚠ Warnings

- All reagents are intended for in vitro "research".
- Reagents have passed Quality Control checks and are expected to react as described above.
 - Should controls samples yield unexpected results, reagents are not to be considered reliable.
- Do not mix reagents from different lot numbers.
- Calibrators dilutions are stable for a few hours.

Technical Information

Please call Customer Support Service or visit www.newscientific.com.

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BIOMEDICINA – PATOLOGIA CLINICA - INFORMATICA

C1Q - Assay Parameters set up on BNA - BN100

Reaction Volumes		Measurement Parameters		Calibration	
Sample	50 µl	Measure Type	Fixed Time	Set Points (Levels)	6
Buffer (*)	240 µl	First Measure	6 sec	First dilution	1:5
Reagent	40 µl	Last Measure	6 min	Sample dilution	1:5
(*) = N reaction buffer					

C1Q - Assay Parameters set up on BNII - ProSpec

Reaction Volumes			Measurement Parameters			Calibration	
Sample	pre reaction	reaction	Last Measure	pre reaction	reaction	Set Points (Levels)	6
	no	µl		50	µl		
Buffer (*)	no	µl	Measure Type	Fixed Time		First dilution	1:5
	no	µl		240	µl		
Reagent	no	µl				Sample Dilution	1:5
	no	µl		40	µl		
(*) = N reaction buffer							
Note Complete BNII assay parametres are on document code: NSC.C1Q.AP.BNII Valid revision A411							

*ImmunoNefelometria su Analizzatori serie BN™
ImmunoNephelometry on "BN™ series" Analyzers
ImmunoNephelometrie auf Analysegeräten Serie BN™*


REF

K.BNA.C1Q
KR.BNA.C1Q

IVD

Reagenti per uso professionale
Reagents for professional use
Reagenzien für Berufsgebrauch

C1Q – Calibration Curve Example on BNII

BNII 2.30

Gio, 4 nov 2004, 15:17

Page 1 of 1


Please note: the value of the curves are only indicative

*Complemento C1q
Complement C1q
Complement C1q*

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BIOMICINA - PATOLOGIA CLINICA - INFORMATICA

C1Q - Assay Parameters set up on BNA - BN100

Reaction Volumes			Measurement Parameters			Calibration	
Sample	50	µl	Measure Type	Fixed Time	Set Points (Levels)	6	
Buffer (*)	240	µl	First Measure	6 sec	First dilution	1:5	
Reagent	40	µl	Last Measure	6 min	Sample dilution	1:5	
(*) = N reaction buffer							

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C1Q - Assay Parameters set up on BNII - ProSpec

Reaction Volumes				Measurement Parameters				Calibration		
Sample	pre reaction	reaction	50	Last Measure	pre reaction	reaction	360	Set Points (Levels)	6	
	no	µl			no	sec				
Buffer (*)	no	µl	240	Measure Type	Fixed Time		First dilution	1:5		
Reagent	no	µl	40						Sample Dilution	1:5
(*) = N reaction buffer										
Note Complete BNII assay parametres are on document code: NSC.C1Q.AP.BNII Valid revision A411										

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New Scientific Company S.r.l. – Via Dante Alighieri, 35 – I 20032 Cormano (MI) – Italy