Valores Asignados	<b>REF</b> W.IGE.STD-K	
	<b>LOT</b> J091	
Assigned Values		
Calibrator Level		IU/ml
Calibrator 1	IGE	0
Calibrator 2	IGE	47
Calibrator 3	IGE	93
Calibrator 4	IGE	186
Calibrator 5	IGE	465
Calibrator 6	IGE	930

BIOMEDICINA - PATOLOGIA CLINICA - INFORMATICA

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**INTENDED USE**

This product line is for use as Calibrators in the determination of:

- **IgE Immunoglobulins (Total IgE) (IGE)**

carried out with NSC's kits.

**METHOD**

The method is based on the Immunoprecipitation reaction in liquid phase (turbidimetry) enhanced with latex particles.

For further information, see the instructions leaflet included in the kits.

**CONTENTS - COMPONENTS**

- |                             |                      |               |
|-----------------------------|----------------------|---------------|
| • <i>IgE - Calibrator 1</i> | <i>(W.IGE.STD-1)</i> | 1 vial - 1 ml |
| • <i>IgE - Calibrator 2</i> | <i>(W.IGE.STD-2)</i> | 1 vial - 1 ml |
| • <i>IgE - Calibrator 3</i> | <i>(W.IGE.STD-3)</i> | 1 vial - 1 ml |
| • <i>IgE - Calibrator 4</i> | <i>(W.IGE.STD-4)</i> | 1 vial - 1 ml |
| • <i>IgE - Calibrator 5</i> | <i>(W.IGE.STD-5)</i> | 1 vial - 1 ml |
| • <i>IgE - Calibrator 6</i> | <i>(W.IGE.STD-6)</i> | 1 vial - 1 ml |

**Form:** liquid, ready for use.

**Stability:** 18 months at +2°C ... +8°C.

**INFORMATION - PREPARATION**

**Origin:** Human serum.

**Concentration**

See table on the back.

Values have been assigned by "Internal Standardisation", using calibrators based on international standards (WHO 1981).

**PRECAUTIONS AND HAZARDS**

- **Preservative:** *Sodium Azide*.

- **Handling:** Sodium Azide is toxic.

All the reagents contain Sodium Azide <0,1% (w/v). Even if Sodium Azide is not dangerous at this concentration, make sure to take the necessary precautions to prevent accidental swallowing.

- **Waste:** Sodium Azide can react with lead and copper plumbing to form highly explosive metal azides.

- **Calibrators and Controls – Serum and Plasma of human origin.**

- **Calibrators and Controls Supplier statement:**

Serum or plasma of human origin has been tested by an "FDA – USA" approved method and found negative for the presence of Hepatitis B surface antigen (HBsAg), HCV antibodies, and HIV-1, 2 antibodies.

- **Handling and waste:**

Because no test method can completely exclude the presence of infectious agents, all materials of human origin should be handled as recommended for any potentially infectious human serum or blood specimen following the precaution recommended for biohazard materials.

**WARNINGS**

- All reagents are intended only for *In Vitro* tests.
- Reagents have passed quality control checks and are expected to react as described in this Instructions. If controls do not yield the predicted reaction, all the reagents would not be considered reliable.
- Do not mix reagents from different lot kits.
- Calibrator dilutions are stable for a few hours.

**TECHNICAL INFORMATION**

Contact Customer Support Service or visit [www.newscientific.com](http://www.newscientific.com).