Assigned Values	REF W.SPC.IS1	
	LOT F051	
ANALYTE		mg / dl
IgA Immunoglobulins	IGA	222
IgG Immunoglobulins	IGG	1 160
IgM Immunoglobulins	IGM	96
Kappa Light Chains Kappa Immunoglobulins	KAP	277 955
Lambda Light Chains Lambda Immunoglobulins	LAM	152 523

BIOMEDICINA - PATOLOGIA CLINICA - INFORMATICA

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INTENDED USE

This product line is for use as Calibrators in the determination of:

- IgA Immunoglobulins (IGA)
- IgG Immunoglobulins (IGG)
- IgM Immunoglobulins (IGM)
- Kappa Light Chains (B&F) - Kappa Immunoglobulins (KAP)
- Lambda Light Chains (B&F) - Lambda Immunoglobulins (LAM)

carried out with NSC's kits.

METHOD

The method is based on the Immunoprecipitation reaction in liquid phase (turbidimetry) with specific antibodies.

For further information, see the instructions leaflet included in the kits.

CONTENTS - COMPONENTS

- *Immunoglobulins - Calibrator* (W.SPC.IS1) 1 vial of 1 ml

Form: liquid, ready for use.

Stability: 12 months at +2°C ... +8°C.

INFORMATION - PREPARATION

Origin: Human serum.

IGA, IGG and IGM values have been assigned based on the *Reference Preparation for Proteins in Human Serum (RPPHS) (*)*.

(*) *RPPHS (Reference Preparation for Proteins in Human Serum)* is a certified reference material developed

jointly through the *College of American Pathologists (CAP)*, the *Bureau Communautaire de Reference* and the *International Federation of Clinical Chemists (IFCC)*.

For further information regarding this reference material please refer to the *RPPHS Lot #91/06 19 (CAP designation)* and *CRM 470 (BCR designation)* product insert.

KAP and LAM values have been assigned based on the *M.M. Lievens* formule (*J.Clin.Chem.Clin.Biochem 1989;27:519-23*).

Concentration

See table on the back.

PRECAUTIONS AND HAZARDS

- Preservative: ***Sodium Azide***.
 - **Handling:** Sodium Azide is toxic.
All the reagents contain Sodium Azide <0,1% (w/v). Even if Sodium Azide is not dangerous at this concentration, make sure to take the necessary precautions to prevent accidental swallowing.
 - **Waste:** Sodium Azide can react with lead and copper plumbing to form highly explosive metal azides.
- Calibrators and Controls – Serum and Plasma of human origin.
 - **Calibrators and Controls Supplier statement:**
Serum or plasma of human origin has been tested by an "FDA – USA" approved method and found negative for the presence of Hepatitis B surface antigen (HBsAg), HCV antibodies, and HIV-1, 2 antibodies.
 - **Handling and waste:**
Because no test method can completely exclude the presence of infectious agents, all materials of human origin should be handled as recommended for any potentially infectious human serum or blood specimen following the precaution recommended for biohazard materials.

WARNINGS

- All reagents are intended only for *In Vitro* assays.
- Reagents have passed quality control checks and are expected to react as described in this Instructions. If controls do not yield the preveded reaction, all the reagents would not be considered reliable.
- Do not mix reagents from different lot kits.
- Calibrator dilutions are stable for a few hours.

TECHNICAL INFORMATION

Contact Customer Support Service or visit www.newscientific.com.