## Italiano

### Scopo

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- Control Complement System Proteins 1 – Medium Value (Y.SPR.CS1.MV) 1 fiala – 1 ml
- Control Complement System Proteins 1 – Low Value (Y.SPR.CS1.LV) 1 fiala – 1 ml

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## English

### Use

This product line is for use in Quality Control in the determination of:

- Complement Component 1q (C1Q)
- C1 (esterase) inhibitor (CEI)
- Properdin factor B (GBG)

carried out with NSC's kits.

### Method

The method used by the kits is based on the reaction of ImmunoPrecipitation in liquid phase with specific antisera. For further information, see the instruction leaflet included in the kits.

### Pack – Form – Stability

- Control Complement System Proteins 1 – Medium Value (Y.SPR.CS1.MV) 1 vial – 1 ml
- Control Complement System Proteins 1 – Low Value (Y.SPR.CS1.LV) 1 vial – 1 ml

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

### Information and Preparation

Human serum, delipidated, defibrinated and filtered to 0,22 µm.

Due to lack of provisions from the Scientific Committees, our "Internal Standardisation" has been applied.

#### *Assigned values*

See table

### ⚠ Precautions and Hazards

Preservative: *Sodium Azide*.

- Caution: Sodium Azide is toxic.

All reagents contain Sodium Azide < 0,1% (w/v). **Sodium Azide is not dangerous at this concentration**.

Careful handling is recommended to avoid accidental swallowing.

- Waste: Sodium Azide may form explosive compounds in metal drainage.

**Calibrators and Controls – Serum and Plasma of human origin**

- Calibrators and Controls Supplier statement Serum or plasma of human origin has been tested by an "FDA – USA" approved method and found negative for the presence of Hepatitis B surface antigen (HBsAg), HCV antibodies, and HIV-1, 2 antibodies.
- Handling and waste Because no test method can completely exclude the presence of infectious agents, all Calibrators and Controls should be handled as recommended for any potentially infectious human serum or blood specimen following the precaution recommended for biohazard materials.

### ⚠ Warnings

- All reagents are intended for *in vitro research*.
- Reagents have passed Quality Control checks and are expected to react as described above. Should controls samples yield unexpected results, reagents and Controls are not to be considered reliable.
- Do not mix reagents from different lot numbers.

### Technical Information

Please call Customer Support Service or visit www.newscientific.com.

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New Scientific Company S.r.l. – Via Dante Alighieri, 35 – I 20032 Cormano (MI) – Italy