## Italiano

### Scopo

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## English

### Use

This product line is for use in Quality Control in the determination of Free Light Chains carried out with NSC's Free Light Chain kits.

### Method

The method used by the Free Light Chain kits is based on the reaction of ImmunoPrecipitation in liquid phase with specific adsorbed antisera for the "hidden" epitopes of the Light Chains of Immunoglobulins.

For further information, see the instruction leaflet included in the Free Light Chains kits.

### List – Pack – Contents – Form – Stability

- |                          |   |             |               |
|--------------------------|---|-------------|---------------|
| <input type="checkbox"/> | <b>Control Free Light Chains kappa</b>  | (Y.UPR.FRK) | 1 vial – 1 ml |
| <input type="checkbox"/> | <b>Control Free Light Chains lambda</b> | (Y.UPR.FRL) | 1 vial – 1 ml |

Form: liquid, ready to use.

Stability: over 12 months at +2°C ... +8°C.

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Due to lack of previsions from the Scientific Committees our "Internal Standardisation" has been applied.

#### *Controls Preparation*

For each type of light chain, samples of urine are collected from patients having Micromolecular Myeloma, selecting those which, by electrophoresis, show a band made up of Bence Jones proteins and only traces of albumin and other proteins. In a sample of this type the Total Proteins quantity can be assumed to be equal to the quantity of Free Light Chains. Total Proteins are measured by the Bradford method.

The sample is centrifuged, dialysed and ultra-filtered. The final product contains Free Light Chains and traces of other proteins in PBS with Sodium Azide < 0,1 % (w/v).

#### *Assigned Values*

See Table.

### ⚠ Precautions and Hazards

Preservative: *Sodium Azide*.

- Caution: Sodium Azide is toxic.
  - All reagents contain Sodium Azide < 0,1% (w/v). *Sodium Azide is not dangerous at this concentration*. Careful handling is recommended to avoid accidental swallowing.
- Waste: Sodium Azide may form explosive compounds in metal drainage.

### ⚠ Warnings

- All reagents are intended for *in vitro research*.
- All kits Reagents, Calibrators and Controls have passed Quality Control checks and are expected to react as described in Kits insert.
  - Should controls samples yield unexpected results, reagents are not to be considered reliable.
- Do not mix reagents from different lot numbers.
- Calibrators dilutions are stable for a few hours.

### Technical Information and References

Please call Customer Support Service or visit [www.newscientific.com](http://www.newscientific.com).

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New Scientific Company S.r.l. – Via Dante Alighieri, 35 – I 20032 Cormanò (MI) – Italy